

Operating Manual **pfm Rotary 3006 EM**

Please read this operating manual carefully before using the instrument for the first time and make certain that you are familiar with the method of operation as well as the instrument's operating elements to ensure that it is used in the correct way.

www.pfmmedical.com



Note

Dear customer,

Please read this operating manual carefully before using the instrument for the first time and make certain that you are familiar with the method of operation as well as the instrument's controls to ensure that it is used correctly.

All the information in this operating manual applies to the pfm Rotary 3006 EM.

Serial number: _____

Please enter the serial number of your instrument here. It can be found on the nameplate.

Part 1 – Safety instructions	5
1.1 Key to symbols	5
1.2 Intended purpose	7
1.3 User group/user profile	7
1.4 Microtome safety devices	8
1.5 Operating the instrument	9
1.6 General precautions	10
Part 2 – Introduction	11
2.1 Views of the pfm Rotary 3006 EM	11
2.2 Technical data for the pfm Rotary 3006 EM	12
2.3 Technical data control panel	12
2.4 Scope of delivery	13
Part 3 – Operating instructions	14
3.1 Transportation, unpacking and setting up	14
3.2 Assembling the accessories	15
3.3 Initial operation	20
3.4 Switching on the instrument	20
3.5 Operating controls	20
3.6 Operating display on the instrument/control panel	21
3.7 Control panel	21
3.8 Menu settings	22
3.9 Selecting Trimming and Coarse Feed modes	22
3.10 Function and adjustment of the Trimming and Coarse Feed modes	22
3.11 Setting the section thickness	23
3.12 Retraction of sample	23
3.13 End position	23
3.14 Sectioning speed	23
3.15 Setting the sectioning speed	23
3.16 Sectioning modes	24
3.17 Starting/stopping the sectioning process (motor-driven)	24
3.18 Sectioning the sample manually/motor-driven	25
3.19 Changing samples/interruption of work	25
3.20 Memo	25
3.21 Remaining feed indication	26
3.22 Section counter/addition	26
3.23 Mounting the object	26
3.24 Standard object clamp	26
3.25 Universal cassette clamp	27
3.26 Object orientation	27
3.27 Adjusting the orientable adapter	27
3.28 Setting the clearance angle	28
3.29 Connecting the foot switch (optional accessories)	28
3.30 Removing/attaching the object mounting device	28
3.31 Handwheel lock	29
3.32 Knife holder	29
3.33 Knife holder base	29

3.34 Knife holder mounting device	30
3.35 Disposable blade holder upper jaw low profile (narrow blade)	30
3.36 Section waste tray	31
3.37 Malfunctions and their correction	31
3.38 Change history	31

Part 4 – Working with the microtome	32
--	-----------

4.1 Preparing, embedding and sectioning the object	32
4.2 Information about how to avoid problems	32
4.3 Trouble-shooting	33

Part 5 – Microtome maintenance, cleaning and care	34
--	-----------

5.1 Maintenance	34
5.2 Cleaning and care	34
5.3 Cleaning and care touch screen/display	34
5.4 Returning the microtome for repair or maintenance	35
5.5 Decontamination	35
5.6 Instrument fuse	35
5.7 Warranty	36
5.8 Declaration of Conformity	36

Part 6 – Disposal of the microtome when put out of use	37
---	-----------

Part 7 – Accessories	38
-----------------------------	-----------



Part 1 – Safety instructions

Please read this operating manual carefully before working with the pfm Rotary 3006 EM. Make certain that you are familiar with the instrument's controls and the method of operation. Please observe the following precautions during each operation of the instrument. Failure to observe these precautions is contrary to the recognised standards of good practice and the intended use of the instrument.

In case of infringement, pfm medical gmbh does not accept any form of liability.

1.1 Overview of symbols and their meanings



Caution! This symbol indicates a possible danger to life and health. Failing to observe these instructions can result in severe health impacts, including life-threatening injuries.

NOTE

Indicates useful information

WARNING

Indicates a dangerous situation



Manufacturer



Date of manufacture



Consult operating manual

NOTE The current electronic operating manual can be found on our website:
<http://www.pfmmedical.com/downloads>



In vitro diagnostic (IVD) medical device



Catalogue number



Serial number



CE marking



Biohazard warning




The instrument must not be disposed of in normal domestic waste.
It must be collected separately.




Keep dry




Fragile, handle with care



Permissible temperature range at the place of use
+10°C - +40°C



Permissible temperature range for storage and transportation
-10°C - +50°C



Permissible humidity range for storage and transportation, and at the place of use
10% - 80% (non-condensing)

1.2 Intended purpose

The pfm Rotary 3006 EM is a modern, fully electronic rotary microtome designed to generate thin sections of human and animal tissue samples of varying hardness, provided the tissue is suitable for sectioning. The instrument can be used for all applications in industry, research and medicine as a diagnostic aid in connection

with in vitro examination as part of the histopathology process. It may only be operated by appropriately trained staff. The sections obtained using the device are intended to assist the pathologist in evaluating the physiological or pathological condition of the tissue samples.

1.3 User group/user profile


- ▶ The pfm Rotary 3006 EM may only be operated by appropriately trained staff who are familiar with the method of operation of the instrument.
- ▶ Before using the instrument for the first time, the operating manual must be read so that the user is familiar with the method of operation and the control elements. Only in this way can proper handling be ensured.

NOTE

The operating manual should be accessible to everyone.

1.4 Microtome safety devices

The instrument is fitted with the following safety devices:



 Finger guard

Finger guard on the knife holder

Before manipulating the knife, knife holder or object, or changing the object, always cover the blade edge with the finger guard.




 Handwheel lock

Handwheel lock

When the object and knife are being mounted, the handwheel lock should be activated to prevent injury during the mounting process. The **M-STOP** field lights up when the handwheel is locked.

 M-STOP




1 Emergency Stop push button

Emergency Stop function

The **E-STOP** field lights up when the Emergency Stop push button is pressed. The emergency stop function can be deactivated by turning the Emergency Stop push button in the direction indicated by the arrow. This must be followed by a reference run (press Start).

E-STOP




1 Handle

Removable handle

The handle can be removed from the side by pressing the button on the front.

WARNING

Caution! For safety reasons, the handle should always be removed for motor-driven sectioning.



1 Handle button

WARNING

Caution! Safety devices must not be removed from the instrument or changed. The proper functioning of the safety devices must be checked by the user before starting work.

1.5 Operating the instrument

- ▶ Take care when handling knives and microtome blades; the blade edge is very sharp and may cause severe injury.
- ▶ When setting up the sample and the knife/blade holder, the handwheel lock should be activated to prevent injuries. The blade edge should be covered with the finger guard.
- ▶ Select the correct blade holder for the microtome blade (Low Profile/High Profile) or knife being used.
- ▶ When setting up the knife/microtome blade in the blade holder, make sure that the blade edge is in the correct position. The blade edge should always be pointing up.
- ▶ The microtome blade should always be completely covered by the pressure plate of the blade holder and not be protruding.
- ▶ If possible, set up the object first before clamping the knife/blade holder into the knife holder.
- ▶ During the sectioning, do not place your fingers between the sample and the knife/microtome blade.
- ▶ A knife/blade holder should always be securely fixed to the microtome so that it cannot fall off.
- ▶ Never try to catch a falling knife/blade holder.
- ▶ Take care during sectioning and when removing sections; the blade edge is very sharp and may cause severe injury.
- ▶ When sectioning, always wear the necessary personal protective equipment (PPE) (safety work shoes/protective gloves/cut-resistant gloves/face mask/protective goggles).
- ▶ After completing the work, always remove the knife/blade holder from the instrument. For the blade holder, ensure that the knife/microtome blade has been removed beforehand.
- ▶ A knife/blade holder that is not in use should always be kept in a safe place (e.g. a knife box) for safety reasons. For the blade holder, ensure that the knife/microtome blade is removed beforehand.
- ▶ Observe the warnings.
- ▶ Warnings on the instrument must not be removed or obscured.
- ▶ Only use suitable accessories and connecting cables recommended/supplied by pfmmmedical.
- ▶ If the instrument fuse needs to be replaced, the parameters specified by the factory must be observed (T 2.0A).
- ▶ Before starting the reference run always move the blade holder to the front position. Otherwise the blade holder may be damaged, because, depending on the position of the object holder, the latter moves downward during the reference run.
- ▶ If unusual noises develop, please contact your responsible pfm service technician or contact pfm medical gmbh directly.
- ▶ The instrument may only be opened and repaired by service technicians trained and authorised by pfmmmedical.
- ▶ No corrosive, acidic or alkaline solutions should be used in the vicinity of the instrument. These could cause corrosion on the instrument and accessories, and damage them.
- ▶ To prevent slip hazards, laboratory shoes should be worn and any paraffin residue should be removed from the microtome and surrounding area.

1.6 General safety precautions

Please take special note of the following warnings:

Beware of explosive gases

The instrument may not be used if dangerous gas concentrations could occur in the surroundings.

Risk of infection

When working with infectious materials, appropriate safety measures must be implemented. The samples used may potentially be infectious. The general laboratory rules regarding protection against infections must therefore be observed.

Risk of radioactivity

When working with radioactive materials, the radiation protection regulations must be followed. Protective clothing must be worn in accordance with the rules and regulations that apply in the respective laboratory for handling radioactively contaminated materials. For the disposal of radioactive waste, the applicable regulations will apply in each instance.

Risk of chemical substances

For better cutting properties, some of the component alloys of this instrument contain small amounts of lead, which is firmly bound in the metal matrix. When used as intended, the lead in the alloys is not released and does not represent any risk to the user. As an additional precaution, we recommend that you wear gloves when changing the accessory parts. You can find additional information in the Declaration on the EU Regulation 1907/2006 (REACH) at: <http://www.pfmmedical.com/downolads>

Avoiding malfunctions

To avoid malfunctions the instrument may not be used close to radio transmitters, e.g. mobile phones. If the instrument malfunctions, please switch it off and contact your distributor.

Caution

Before switching the instrument on for the first time, please check at the installation location that the supply voltage/frequency match the information given on the nameplate.

Protective earthing of the instrument

In order to rule out the risk of electrocution, the instrument must be connected to the earth conductor. The instrument is equipped with a mains power plug with earth conductor. It is very important to check that the mains power outlet used is equipped with an earth conductor, and that it complies with the regulations of the International Electrotechnical Commission (IEC).

Mains voltage precautions

Never remove housing components. The replacement of components and alignment must only be carried out by trained service technicians.

Interference due to electromagnetic interactions


The electromagnetic environment should be assessed before the microtome is operated. Do not use this instrument near strong sources of electromagnetic radiation (for example, unshielded, intentionally operated high-frequency sources) because these can interfere with the proper operation of the instrument. The pfm Rotary 3006 EM fulfils the requirements of the standard IEC 61326-2-6 on immunity and emissions.

WARNING


Caution! Switch off the instrument and pull out the mains plug before removing any part of the housing.

Part 2 – Introduction


2.1 Views of the pfm Rotary 3006 EM



- 1 Display
- 2 Control wheel Settings
- 3 Control wheel Coarse Feed/Trimming
- 4 Control wheel Sectioning Speed
- 5 Emergency Stop push button
- 6 Clamping lever for clearance angle adjustment
- 7 Clamping levers for knife holder
- 8 Clamping lever for object orientation
- 9 Handle button
- 10 Section waste tray



- 11 Finger guard
- 12 Pressure plate
- 13 Knife holder
- 14 Handwheel
- 15 Adjustment screws for object orientation with tangible zero-point
- 16 Handwheel lock
- 17 Clamping lever for pressure plate
- 18 Knife angle scale



2.2 Technical data for the pfm Rotary 3006 EM

- ▶ **Section thickness:** 0.5 - 100 μm
- ▶ **Setting range:** 0.5 - 10 μm in 0.5 μm steps / 10 - 20 μm in 1 μm steps / 20 - 50 μm in 5 μm steps / 50 - 100 μm in 10 μm steps
- ▶ **Trimming mode:** 1 - 500 μm
- ▶ **Coarse feed speed:** 1.0 mm/sec.
- ▶ **Sectioning modes:** manual, continuous, single and step
- ▶ **Sectioning speed max.:** 300 mm/sec.
- ▶ **Retraction:** 0 - 100 μm in 5 μm steps (in manual operation)
- ▶ **Remaining feed indication 1 mm horizontal:** visual and audible
- ▶ **Horizontal range:** 30 mm
- ▶ **Vertical range:** 70 mm
- ▶ **Object orientation:** x- and y-axes universal 12° z-axis up to 360°
- ▶ **Maximum sample size:**
 - ▶ Universal cassette clamp, standard embedding cassettes
 - ▶ Standard object clamp large, sample size 45 x 60 mm
 - ▶ Universal cassette clamp large, super mega embedding cassettes
- ▶ **Power supply:** 100/120/230/240 V AC, 50 - 60 Hz
- ▶ **Operating environment altitude:** up to 2,000 m above sea level
- ▶ **Pollution degree (according to IEC 61010-1):** II
- ▶ **Overvoltage category (according to IEC 61010-1):** II
- ▶ **Fluctuations in mains voltage:** $\pm 10\%$ of the mains voltage
- ▶ **IP protection class:** IP20
- ▶ **Dimensions (W/D/H):** 450 x 530 x 275 mm
- ▶ **Weight:** 38 kg

2.3 Technical data control panel

- ▶ **Dimensions (W/D/H):** 135 x 170 x 55 mm
- ▶ **Height opened:** 150 mm
- ▶ **Weight:** 1.9 kg

2.4 Scope of delivery

pfm Rotary 3006 EM

REF: 0300610

Scope of delivery

Base device

1 section waste tray

1 mains cable

1 protective cover

Without blade holder base and disposable blade holder upper jaw

Without object clamp and adapter

Without foot switch

Without control panel

NOTE

The current electronic operating manual can be found on our website: <http://www.pfmmedical.com/downloads>

Part 3 – Operating instructions

3.1 Transportation, unpacking and setting up the instrument

- ▶ Before unpacking the delivery, please check that it is complete and undamaged.
- ▶ Disassemble the packaging and remove the instrument and accessories from the packaging. The section waste tray is included in the instrument packaging, but the other accessories are packaged separately.
- ▶ Take care when unpacking the instrument. The instrument is very heavy and should be removed from the packaging and transported to the installation location by at least two people for safety reasons.
- ▶ When lifting the instrument out of the packaging and during transportation, the instrument should only be handled at the base plate. Do not carry or lift the instrument at the handwheel or object mounting device.
- ▶ Please read the operating manual carefully before using the instrument for the first time.
- ▶ When setting up or moving the instrument, make sure not to place your fingers between the instrument and the set-up area.
- ▶ The instrument should be positioned on a stable, vibration-free table.
- ▶ The instrument must only be positioned on non-flammable, flat surfaces.
- ▶ Position the instrument so that the handwheel is freely and easily accessible.
- ▶ If other instrument in the vicinity of the microtome generates vibrations, these can interfere with the sectioning process.
- ▶ Position the instrument so that the power plug is always accessible.
- ▶ Ambient temperature at the installation location of the instrument between + 10°C and + 40°C.
- ▶ Relative humidity at the installation location of the instrument max. 80% (non-condensing).
- ▶ Ambient temperature during the transportation and storage of the instrument: between -10° and +50°C.
- ▶ Relative humidity during the transportation and storage of the instrument: max. 80% (non-condensing).

The operator is responsible for checking that the instrument is working properly and safely.


The operator shall assume responsibility for any use of special instrument and materials for operating the instrument.

NOTE

If possible, retain the original packaging for the entire lifetime of the instrument so that, if it requires servicing, you can package the instrument correctly again.

3.2 Assembling the accessories

Assembling the object clamps with/without adapter



- 1 Clamping lever for object orientation
- 2 Adjusting screws for object orientation
- 3 Adapter orientable



Adapter orientable



Universal cassette clamp, orientable




Universal cassette clamp, non-orientable (without adapter)

- ▶ Release and remove the clamping lever for object orientation.
- ▶ Turn the two adjustment screws for object orientation to the tangible zero-point (marked in red).
- ▶ Insert the object clamp with the adapter, paying attention to the position of the hole so that the clamping lever can later be inserted without difficulty.
- ▶ Clamp the object clamp using the clamping lever.
- ▶ After each change, check that the clamping lever is sufficiently clamped.

NOTE

Some object clamps are optionally also available as a non-orientable version. For use on the instrument, the object clamp must be attached to the adapter fixed or adapter orientable. It is then fitted to the instrument as described.


Knife holder




- 1 Finger guard
- 2 Disposable blade holder upper jaw
- 3 Knife holder base

The knife holder consists of the knife holder base and the disposable blade holder upper jaw. The finger guard is fitted to the disposable blade holder upper jaw.


Assembling the knife holder base



- 1 Knife holder base
- 2 Knife holder clamping levers





- 1 Guide rails
- 2 Knife holder mounting device



- 1 Knife holder base

- ▶ Release the knife holder clamping levers.
- ▶ Slide the knife holder base onto the base plate via the two guides and fix it above the knife holder mounting device.
- ▶ Clamp the knife holder mounting device using the knife holder clamping levers.
- ▶ After each change, check that the clamping levers are sufficiently clamped.

Assembling the disposable blade holder upper jaw



- 1 Pressure plate
- 2 Clearance angle adjustment clamping lever
- 3 Pressure plate clamping lever
- 4 Knife angle scale

- ▶ Release and completely remove the clearance angle adjustment clamping lever.
- ▶ Insert the adapter on the disposable blade holder upper jaw into the knife holder base. Take note of the position of the hole, so that the clamping lever can later be re-inserted without difficulty.
- ▶ Set the clearance angle to 10° using the knife angle scale.
- ▶ Clamp the clearance angle adjustment clamping lever.
- ▶ After each change, check that the clamping lever is sufficiently clamped.

Connecting the foot switch



1 Foot switch connection

Foot switch

- ▶ The colour-marked (blue) connection for the foot switch is located on the back of the instrument.
- ▶ Connect the foot switch cable to the foot switch socket.

WARNING

Caution! Please observe the guide on the foot switch socket and do not turn the plug in the connection socket.

Connecting the control panel



1 Control panel connection



Control panel

- ▶ The colour-marked (green) connection for the control panel is located on the back of the instrument.
- ▶ Connect the control panel cable to the control panel socket.

WARNING



Caution! Please observe the guide on the foot switch socket and do not turn the plug in the connection socket.

Section waste tray



The section waste tray is pushed onto the microtome via two guides on the inside and then positioned so that it is felt to latch.

Fitting the microtome blade




Feather® microtome blade S35

- ❶ Finger guard
- ❷ Pressure plate
- ❸ Blade support
- ❹ Pressure plate clamping lever

- › Fold down the finger guard.
- › Release the clamping lever for the pressure plate.
- › Remove the microtome blade from the dispenser and carefully slide it under the pressure plate from the side.
- › It is easier to introduce the microtome blade by gently pressing on the pressure plate. This forms a small gap between the blade support and the pressure plate, making it easier to introduce and position the microtome blade.
- › Make sure that the microtome blade is positioned parallel to the pressure plate.
- › Clamp the pressure plate using the clamping lever.
- › After each change, check that the clamping lever is sufficiently clamped.


Fitting a standard knife



- ❶ Finger guard
- ❷ Screws
- ❸ Knurled screw

- › Fold down the finger guard.
- › Loosen the two screws on the standard knife holder.
- › Insert the knife into the knife holder from the side.
- › The height of the knife can be adjusted and it can be aligned parallel to the finger guard using the two knurled screws at the side.
- › Clamp the knife using the two screws.
- › After each change, check that the knife is sufficiently clamped.

3.3 Initial operation



- 1 Connection socket
- 2 instrument fuse
- 3 Power switch

- ▶ Before switching on the microtome for the first time, you must make sure that the supply voltage and frequency specified on the nameplate match the mains voltage specifications at the installation location.
- ▶ Only connect the instrument to earthed wall sockets.
- ▶ Connect the mains power cable to the connection socket at the back of the instrument.
- ▶ The insert for the instrument fuse (T2.0A) is situated next to the connection socket.
- ▶ Now switch on the instrument using the power switch.
- ▶ 1=ON, 0=OFF

WARNING

Caution! Only use the supplied mains power cable with a ferrite core.

Caution! Do not operate the instrument if the mains power cable is damaged.

Caution! Do not use an extension cable for the connection.

Caution! If the instrument is not used with the supplied mains power cable, the user must carry out a safety test with the new mains power cable before the instrument is used for the first time.

3.4 Switching on the instrument



After switching on the instrument an audible signal can be heard, and "REFERENCING" is shown on the display. The reference run is started by pressing **START**. After this, the display automatically goes back to the control panel display.

WARNING

Caution! Before starting the reference run always move the knife holder to the front position. Otherwise the knife holder may be damaged, because, depending on the position of the object holder, the latter moves downward during the reference run.

Caution! Do not operate the instrument if the display has no function, shows illogical or incorrect values or is damaged.

3.5 Operating controls

The pfm Rotary 3006 EM is equipped with a modern touch screen system. The control panels and displays are clear and arranged by function. The user can operate the instrument either via the operating display on the instrument or via a separate control panel.

This allows for simple and comfortable operation. When the control panel is connected, the instrument can only be controlled from the control panel. The control panel is not included in the general scope of delivery and must be ordered separately.

3.6 Operating display on the instrument/control panel



- 1 Emergency stop
- 2 Modes
- 3 Section thickness
- 4 Trimming
- 5 Trimming automatically/manually
- 6 Retraction
- 7 Section counter
- 8 Remaining feed
- 9 Speed
- 10 Cutting Window on/off
- 11 Menu
- 12 Start
- 13 Delete
- 14 Memo
- 15 M-STOP
- 16 Control wheel Settings
- 17 Control wheel Coarse Feed/Trimming
- 18 Control wheel Sectioning Speed

3.7 Control panel



If the pfm Rotary 3006 EM is to be operated via the separate control panel, it must be connected to the instrument before switching it on. The connection for the control panel is located on the back of the instrument.

3.8 Menu settings



Press the **MENU** button to get to a submenu with the following subfunctions:

Buzzer

Press the **BUZZER** button to switch the buzzer on or off.

Language

Press the **LANGUAGE** button to select between different languages.

Brightness

Press the **BRIGHTNESS** button to adjust the brightness of the display. There are three levels in total.

Reset/factory default settings

Press the **RESET** button to reset the instrument to factory defaults.

3.9 Selecting Trimming and Coarse Feed modes

By pressing the Trimming/Coarse Feed control wheel you can switch between Trimming and Coarse Feed.

3.10 Function and adjustment of the Trimming and Coarse Feed modes

By pressing the Trimming/Coarse Feed control wheel you can switch between Trimming and Coarse Feed. Press the **TRIMMING** button to select the required value by turning the Settings control wheel. Press the **TRIMMING** button again or press the Settings control wheel to save the set value.

Trimming function

The Trimming function can be used to select either manual or automatic sectioning of the sample. The required function is activated by pressing the **MAN** or **AUT** button. When the **TRIMMING** **AUT** function has been activated, the **AUT** button blinks.

Adjusting the trimming thickness during the trimming

It is also possible to change the preset trimming thickness during the trimming by turning the Settings control wheel.

Automatic trimming (Trimming AUT)

In automatic trimming the function is activated by briefly turning the Trimming/Coarse Feed control wheel forwards and deactivated again by briefly turning it forwards again. If the function is deactivated the infeed automatically reverts to the preselected section thickness.

WARNING

Caution! The advance via **Trimming AUT** is interrupted only when the function is deactivated again using the control wheel (turning it forwards).

Caution! Always adjust the trimming thickness setting to the size of the sample.

Manual trimming (Trimming MAN)

In manual trimming, infeed occurs only when the Trimming/Coarse Feed control wheel is turned forwards during sectioning. If the control wheel is released, the infeed automatically reverts to the preselected section thickness.

Coarse feed functions

The coarse feed function is intended to quickly change the distance from the sample to the knife or away from the knife. In the Coarse Feed mode the advance feed only works when the Trimming/Coarse Feed control wheel is turned forward.

3.11 Setting the section thickness

Press the **THICKNESS** button to select the required section thickness by turning the Settings control wheel. Press the **THICKNESS** button again or press the Settings control wheel to save the set value.

Adjusting the section thickness during sectioning

It is also possible to change the preset section thickness during the sectioning by turning the Settings control wheel.

3.12 Retraction of sample

The retraction function is activated by pressing the **RETRACTION** button. Turn the Settings control wheel to select the required value. Press the **RETRACTION** button on the display again or press the Settings control wheel to save

the set value. When the value is 0, retraction is deactivated. Retraction is activated during the upwards motion in order to protect the knife/blade and object from damage.

3.13 End position

If you stop the sectioning process, the instrument remains in the end position (Continuous mode “Up”, Single mode “Up”). For stopping the sectioning operation, see 3.17.

3.14 Sectioning speed

The sectioning speed can be adjusted in the area of the cutting window. To activate this function, press the Cutting Window button. This enables you to reduce the sectioning speed during the sectioning process.

3.15 Setting the sectioning speed

By turning the Sectioning Speed control wheel, the sectioning speed can be continuously adjusted. Turning it clockwise makes sectioning faster, and turning it anticlockwise makes sectioning slower.

WARNING

CAUTION! Always adjust the sectioning speed to the composition of the sample being sectioned. The harder the sample, the slower the sectioning speed.

3.16 Sectioning modes (manual/motor-driven)

The instrument can be operated either manually or motor-driven. You have one manual mode and three different motor-driven modes at your disposal.

Selecting the mode

The desired operating mode can be selected by directly pressing the icons on the display.

Modes

Manual

The sectioning movement is controlled manually by turning the handwheel. The rotary motion is clockwise. The section is made as the object moves downwards. In manual sectioning mode **MAN** the instrument can be switched to motor-driven **CON** continuous operation using the start button, the Sectioning Speed control wheel or the foot switch. If the instrument is stopped, the operating mode automatically changes back to the manual sectioning mode **MAN**.

(This function is only available from version 4.4.3 onwards. For older versions before 4.4.3 please contact pfm medical gmbh.)

Motor-driven CONTINUOUS

Continuous sectioning operation until the instrument is stopped. After that, the object stops in the selected end position (Up).

Motor-driven SINGLE

Single stroke After that, the object stops in the selected end position (Up).

Motor-driven STEP

Continuous operation by permanently pressing the **START** button, the Sectioning Speed control wheel, or by permanently using the foot switch. On release, the object stops immediately.

WARNING

Caution! In manual sectioning mode **MAN** using the start button, the Sectioning Speed control wheel or the foot switch changes the instrument to motor-driven continuous operation **CON** and the instrument starts immediately.

3.17 Starting/stopping the sectioning operation (motor-driven)

The motor-driven sectioning operation can be started or stopped by pressing the **START** button, by pressing the Sectioning Speed control wheel or by using the foot switch.

3.18 Sectioning the sample manually/motor-driven

Manual

- ▶ Select the **MAN** mode.
- ▶ Release the handwheel lock.
- ▶ Coarsely advance the sample using the Trimming/Coarse Feed control wheel in the Coarse Feed mode.
- ▶ Change to Trimming mode (AUT/MAN) by pressing the Trimming/Coarse Feed control wheel and trim the sample by turning the handwheel.
- ▶ If automatic trimming is deactivated (TRIMMING AUT) or if the Trimming/Coarse Feed control wheel is released (TRIMMING MAN) the infeed automatically reverts to the preselected section thickness.

Motor-driven CONTINUOUS (continuous operation)

- ▶ Select the **CON** mode.
- ▶ Release the handwheel lock.
- ▶ To coarsely advance the object towards the knife, select the Coarse Feed mode using the Trimming/Coarse Feed control wheel.
- ▶ Move the object towards the knife using the Coarse Feed mode.
- ▶ Select the appropriate sectioning speed using the Sectioning Speed control wheel.
- ▶ Start the motor-driven sectioning operation.
- ▶ Change to Trimming mode (AUT/MAN) by pressing the Trimming/Coarse Feed control wheel and trim the sample.
- ▶ If automatic trimming is deactivated (TRIMMING AUT) or if the Trimming/Coarse Feed control wheel is released (TRIMMING MAN) the infeed automatically reverts to the preselected section thickness.

WARNING

Caution! The handle should be removed for motor-driven operation.

Motor-driven SINGLE (single stroke)

- ▶ Select the **SING** mode.
- ▶ Release the handwheel lock.
- ▶ To coarsely advance the object towards the knife, select the Coarse Feed mode using the Trimming/Coarse Feed control wheel.
- ▶ Move the object towards the knife using the Coarse Feed mode.
- ▶ Select the appropriate sectioning speed using the Sectioning Speed control wheel.
- ▶ Start the motor-driven sectioning operation.
- ▶ The motor is set in motion for one stroke and afterwards stops in the end position.

Motor-driven STEP (step mode)

- ▶ Select the **STEP** mode.
- ▶ Release the handwheel lock.
- ▶ To coarsely advance the object towards the knife, select the Coarse Feed mode using the Trimming/Coarse Feed control wheel.
- ▶ Move the object towards the knife using the Coarse Feed mode.
- ▶ Select the appropriate sectioning speed using the Sectioning Speed control wheel.
- ▶ Start the motor-driven sectioning operation.
- ▶ The motor is set in motion for as long as the **START** button, the Sectioning Speed control wheel, or the foot switch is pressed.
- ▶ On release, the motor stops immediately.

3.19 Changing samples/interruption of work

Changing samples

- ▶ When the instrument is stopped, the motor stops in the end position (See Mode).
- ▶ Use the Trimming/Coarse Feed control wheel to change to the Coarse Feed mode and establish enough distance between the object and the knife.
- ▶ Lock the handwheel.
- ▶ Cover the blade edge with the finger guard.
- ▶ Change the sample.

Interruption of work

- ▶ When the instrument is stopped, the motor stops in the end position. (See Mode)
- ▶ Use the Trimming/Coarse Feed control wheel to change to the Coarse Feed mode and establish enough distance between the object and the knife.
- ▶ Lock the handwheel.
- ▶ Cover the blade edge with the finger guard.

3.20 Memo

By pressing the **MEMO** button (acoustic signal), the sample's present horizontal position is saved.

By pressing the **MEMO** button again, the saved sample position is restored.

3.21 Remaining feed indication



- 1 Rear end position
- 2 Front end position

REMAINING FEED 00000

The remaining feed is indicated acoustically and optically. The potential remaining feed in μm is shown optically on the display. In the reference position "Rear" the value 30000 is shown on the display. In the end position "Front" the value 0 is shown. An acoustic signal sounds between the positions 29000 and 30000 and between 1000 and 0 to indicate that only 1000 μm (1 mm) of feed remains.

WARNING

Caution! A red marker on the object mounting device also indicates that the front or rear end position of the horizontal stroke has been reached.

3.22 Section counter/addition

A section counter is used to add up the sections performed. After each turn of the handwheel, the counter increases by one or by the section thickness in μm .

By pressing the **DEL** button, the section counter/addition is set to zero.

3.23 Mounting the object



- 1 Adapter orientable

The object is mounted using different clamping systems. An orientable adapter allows easy alignment of the object to the sectioning level. The orientable adapter can be screwed to the appropriate clamp if needed.

WARNING

Caution! It is imperative that the handwheel be locked when working in the object orientation area.

3.24 Standard object clamp



The standard object clamp is used for paraffin or plastic blocks of square or rectangular shapes.

WARNING

Caution! The contact surfaces on the cheeks of the object clamps must be free of any residual paraffin. Otherwise proper clamping cannot be guaranteed.

3.25 Universal cassette clamp




The universal cassette clamp, which is operated via the clamping lever, is used for mounting standard embedding cassettes. The clamping lever must be pulled upward to insert or withdraw an embedding cassette.

WARNING

Caution! The contact surfaces on the cheeks of the object clamps must be free of any residual paraffin. Otherwise proper clamping cannot be guaranteed.

3.26 Object orientation




In order to move the object to the required position, the clamping lever must first be released (red marker point up). The clamping lever is set to the clamp position as soon as the required position is reached. The set position is thus fixed.

Zero-point setting

The two adjustment screws must be turned as far as the red markers. For ease of orientation, this position is tangible.

- 1 Clamping lever
- 2 Object orientation adjustment screws

3.27 Adjusting the orientable adapter




To adjust the clamping of the orientable adapter, first insert a screwdriver in the hole in the adapter to fix it firmly. Unscrew the nut with a 13 mm open spanner. Then use an 8 mm open spanner to screw the threaded brass bolt inwards or outwards as required and tighten the nut again with the 13 mm spanner.

To test the clamping, insert the adapter in the guide bore and tighten it with the clamping lever. If the clamping is not adequate, the process must be repeated and the brass bolt adjusted in the required direction. It must be possible to fix the clamping lever at an angle of about 45° in both directions.

- 1 Adapter hole
- 2 Nut
- 3 Threaded brass bolt

3.28 Setting the clearance angle



- 1 Knife angle scale
- 2 Clamping lever for clearance angle adjustment

On the right side of the knife holder there is a scale for adjusting the clearance angle. Release the clamping lever on the left side of the knife holder and set the desired clearance angle, then fix the knife holder in the selected position again using the clamping lever.

3.29 Connecting the foot switch (optional accessory)

The connecting socket for the foot switch is on the back of the instrument. The foot switch can be used to operate the motor-driven sectioning process (Start/Stop). In addition, it has a similar function to the Emergency Stop push button. If the foot switch is pressed firmly, the Emergency Stop function is activated.

WARNING


Caution! Please do not turn the plug in the connection socket. Please observe the guides on the foot switch socket.

WARNING

Caution! In manual sectioning mode **MAN** using the start button, the Sectioning Speed control wheel or the foot switch changes the instrument to motor-driven continuous operation **CON** and the instrument starts immediately.

(This function is only available from version 4.4.3 onwards. For older versions before 4.4.3 please contact pfm medical gmbh.)

3.30 Removing/attaching the object mounting device




- 1 Clamping lever for object orientation
- 2 Adjustment screws for object adapter
- 3 Adapter orientable



To change the object mounting device, release the clamping lever and remove it from the guide bore. Both adjusting screws must be turned as far as the red markers. Then the entire object clamp together with the orientable adapter can be removed upwards. The new object mounting device is fitted in reverse order. Please observe the guide bore position in the adapter to allow reinsertion of the clamping lever. After each change, check that the clamping lever is sufficiently clamped. If it is not properly clamped, vibrations etc. can occur while sectioning.

3.31 Handwheel lock




- 1 Handwheel
- 2 Handwheel lock

The handwheel can be locked within 360°. To lock it, the slide of the handwheel lock is pushed all the way back. To release the lock, the slide must again be pushed (in the opposite direction). The **M-STOP** field lights up when the handwheel is locked.

WARNING

Caution! After completing the sectioning, the handwheel should always be locked using the handwheel lock. Take care that the slide of the handwheel lock is in the correct position to ensure that it is locked correctly.

3.32 Knife holder




- 1 Finger guard
- 2 Disposable blade holder upper jaw
- 3 Knife holder base

The knife holder consists of the knife holder base and the disposable blade holder upper jaw. For better cleaning of the disposable blade holder upper jaw as well as the pressure plate and the contact surfaces, the finger guard on both sides is “only” clipped on and therefore removable.


WARNING

Caution! The finger guard must be used and the handwheel locked during all work on the knife holder. The position of the clamping lever on the knife holder can be adjusted as required (360°).

3.33 Knife holder base



- 1 Knife holder base




- 1 Knife holder base
- 2 Knife holder clamping levers

The knife holder base is used for mounting the disposable blade holder upper jaw and is fixed to the base plate of the microtome above the knife holder mounting device. The knife holder base can be released or clamped with the knife holder clamping levers (on the left and right).

3.34 Knife holder mounting device




- 1 Knife holder mounting device
- 2 Knife holder clamping levers



- 1 Scale


The knife holder base is fixed to the base plate of the microtome above the knife holder mounting device. The knife holder mounting device is operated with the knife holder clamping lever. The mounting device and the precision guide must be free of any residual paraffin and should be cleaned regularly. For optimal positioning of the knife holder to the sample, the two guides are provided with a scale. The knife holder base can be released or clamped using the clamping levers on the left and right.

3.35 Disposable blade holder upper jaw low profile (narrow blade)




- 1 Pressure plate
- 2 Clearance angle adjustment clamping lever
- 3 Pressure plate clamping lever
- 4 Knife angle scale

In order to fix the microtome blade in the disposable blade holder upper jaw with the pressure plate problem-free, the gap between the pressure plate and the mounting rail must first be opened. To do this, move the clamping lever for the pressure plate as far as it will go and insert the microtome blade. The lower end of the pressure plate should be pressed during the insertion of the microtome blade, for optimal use of the gap opening. The microtome blade is also clamped using the clamping lever. The clearance angle is adjusted using the scale on the knife holder base and the index on the disposable blade holder upper jaw. The clearance angle adjustment is released or clamped using the clamping lever. Under normal circumstances, a clearance angle setting of 10° is recommended for the pfm Rotary 3006 EM. The upper jaw of the disposable blade holder is optionally available either as a low profile blade holder for narrow blades or as a high profile blade holder for wide blades.



- 1 Disposable blade holder upper jaw Low Profile (narrow blade)



- 1 Disposable blade holder upper jaw High Profile (wide blade)

3.36 Section waste tray



The section waste tray is anodised and therefore easy to clean. It is pushed onto the microtome via two guides on the inside and then positioned so that it is felt to latch.

3.37 Malfunctions and their correction

Defect	Possible cause	Elimination
No display after switching on the instrument	Mains cable not correctly connected	Check the mains cable connection
	Mains cable defective	Check the mains fuse and replace if necessary
	Instrument fuse defective	Change the instrument fuse
Instrument does not start	Emergency Stop push button has been activated	Deactivate Emergency Stop push button
	Handwheel lock is activated	Release the handwheel lock
No infeed during sectioning or trimming of the sample	The object has reached the front position	Move the object backwards using the coarse feed function
Foot switch does not work	Foot switch not connected correctly	Check connection to foot switch
	Foot switch defective	Check foot switch and replace if necessary
Handwheel is hard to turn	Ball bearing defective	Replace the entire drive claw
The clamping eccentric cannot be fully tightened	Clamping of the adapter is not correctly adjusted	Adjust clamping of adapter

3.38 Change history

State of revision	Change
Rev.2023-03-13	Design reworked
Rev.2024-01-01	Corporate identity guidelines

Part 4 – Working with the microtome

4.1 Preparing, embedding and sectioning the object

The objects must be properly prepared and embedded in suitable media before sectioning. Different forms and rectangular samples as well as the embedding cassette are used for embedding. Depending on the embedding form, the object is mounted in the standard object clamp or in

the universal cassette clamp and then trimmed. Using the coarse feed (Trimming/Coarse Feed control wheel), the object can be coarsely approached to the blade so that there is a first contact between the object and the knife blade. Sectioning can be accelerated using the trimming function.

4.2 Information about how to avoid problems

Sample preparation

When preparing the sample, the correct fixation and embedding medium should be selected, and drainage and infiltration times should be observed.

Sample temperature

Heating of a paraffin sample through direct sunlight or other heat sources should be avoided, otherwise the paraffin will become soft and it will no longer be possible to section the sample.

Clamping screws

Tighten all clamping screws on the knife holder and object holder.

Clearance angle

The knife clearance angle must be selected depending on the sample.

Sectioning speed

The harder the material, the slower the sectioning speed.

Selection of the FEATHER® microtome blades

Information can be obtained directly from pfm medical gmbh.

service@pfmmedical.com

T +49 (0)2236 9641-220

F +49 (0)2236 9641-51

www.pfmmedical.com

4.3 Trouble-shooting

Defect	Possible cause	Elimination
Thick/thin sections	Microtome blade is insufficiently clamped in the knife holder	Check the clamping and clean
	Knife holder is insufficiently clamped	Check the clamping and clean
	Object holder is insufficiently clamped	Check the clamping and clean
	Incorrect knife holder clearance angle	Check, adjust the clearance angle (10°)
Handwheel is hard to turn	Handwheel lock	Check operation/latching
Vibrations while sectioning	Object holder is insufficiently clamped	Check the clamping and clean
	Knife holder/pressure plate is insufficiently clamped	Check the clamping and clean
	Microtome blade in knife holder	Check, adjust position, clean
Section compressed/wrinkled	Incorrect knife holder clearance angle	Check, adjust the clearance angle (10°)
	Block too warm	Cool down the block
Ribbon sectioning is not possible	Incorrect knife holder clearance angle	Check, adjust the clearance angle (10°)
	Knife holder/pressure plate is insufficiently clamped	Check the clamping and clean
Crooked sections	Lateral edges of the block are not parallel	Move them into a parallel position
	Blade edge and block are not parallel	Align the block
	Pressure plate is insufficiently clamped	Check the clamping and clean
Section rolls up	Incorrect knife holder clearance angle	Check, adjust the clearance angle (10°)
	Microtome blade is still too sharp	Longer trimming time required
Section too thick	Object holder is insufficiently clamped	Check the clamping and clean
	Knife holder/pressure plate is insufficiently clamped	Check the clamping and clean
	Block too warm	Cool down the block

Part 5 – Maintenance, cleaning and care

5.1 Maintenance

In order to guarantee the sectioning quality and operation of the rotary microtome, annual routine maintenance should be carried out by a trained service technician authorised by pfmmedical. pfm medical gmbh offers maintenance to ensure that your instrument always

remains in faultless condition. Details can be obtained from your local sales representative or directly from pfm medical gmbh. We recommend an annual electrical safety inspection.

5.2 Cleaning and care

It is recommended that the instrument be cleaned by the user depending on the degree of dirt contamination and frequency of operation. Before cleaning, the instrument must be switched off at the power switch and the mains plug must be disconnected. Mild household cleaning agents

are suitable for the cleaning. Do not use any aggressive cleaning agents or solvents as these could damage any parts made of plastic. No liquids should penetrate the inside of the instrument. The cross roller guides are protected against dirt and do not need to be cleaned.

WARNING

Caution! When handling disinfectants and cleaning agents, always observe the manufacturer's specifications regarding safety instructions, dilution and exposure time.

5.3 Cleaning and care touch screen/display

The display and touch screen are made of plastic and should not be touched with hard, sharp or pointy objects. Do not use any aggressive cleaning agent or solvents for the cleaning. No glass cleaners or detergent should be

used as they contain aggressive soap solutions that attack and damage the surface. We recommend drying the touch screen in as dry a manner as possible; it is best to use a commercially available microfibre cloth.

WARNING

Caution! The display should not be sprayed with a cleaning agent or disinfectant because there is a danger that liquid will get between the films or behind the display, and damage the display.

Caution! When handling disinfectants and cleaning agents, always observe the manufacturer's specifications regarding safety instructions, dilution and exposure time.

5.4 Returning the microtome for repair or maintenance

Generally, repair or maintenance work is carried out on site. If this is not possible for any special reasons, the equipment can be sent to pfm medical gmbh. Before you send a piece of equipment to pfmmedical for a service, it must be cleaned and decontaminated in an appropriate manner. This also applies to equipment undergoing a service on site. Contaminated equipment pose a potential health hazard for anyone coming into contact with them due to infectious agents, pathogens or medicinal products. The legal requirements for the handling of contaminated instruments are based, among others, on the Occupational Safety and Health Act (ArbSchG), the Biological Agents Ordinance (BioStoffV), the Occupational Safety Act (ASiG) and the Accident Prevention Regulations (Unfallverhütungsvorschriften). To make sure that the decontamination was performed, please complete the decontamination certificate carefully. You can find the current certificate on our website: <http://www.pfmmedical.com/downloads>

The service measures can only be carried out if a certificate is fully completed in advance and enclosed with the instrument. If it transpires that the instrument is a potential hazard, we reserve the right to return the instrument to the sender immediately at the sender's expense.

WARNING

Caution! Please only send the instrument in the original packaging together with the mains cable. Please contact pfmmedical if you no longer have the original packaging.


Caution! Knives/blade holders should always be packed in the appropriate box. For the blade holder, ensure that the knife/microtome blade has been removed beforehand.

5.5 Decontamination

For the regular disinfection of the work surfaces and other surfaces of the microtome, we recommend using single-use disinfectant wipes, for example those recommended by

the RKI (Robert Koch Institute) or the VAH (Association for Applied Hygiene) list: "Mikrozid AF" and "Bacillol AF".

5.6 Instrument fuse



- 1 Connection socket
- 2 Instrument fuse
- 3 Power switch

The instrument fuse is located on the back of the instrument, and it is integrated into the power switch. To replace the fuse, the instrument must be switched off at the power switch and the mains plug must be disconnected. Open the small protective cover with a small flat screwdriver and remove the fuse holder. The fuse can then be replaced by a new one.

WARNING

Caution! Before replacing the fuse, the instrument must be switched off and disconnected from the mains power.

5.7 Warranty

We have checked and tested this instrument thoroughly. However, in the event of any defect, apart from minor defects, occurring in a piece of instrument or spare part distributed by us up to one year from delivery, we undertake to remedy the defect to the extent specified below, without prejudice to any statutory or other contractual claims, even where the defect was not present on delivery to the customer, provided that there has been no unauthorised interference with the instrument by third parties. This does not include any assurance of the instrument being free of wear and tear, or any guarantee of the condition or durability of the instrument.

The remedying of defects under this undertaking shall be carried out, at our discretion, by the replacement or repair of defective parts. Spare parts may be supplied for up to 5 years after the delivery of the instrument. Any warranty claims on the part of the customer beyond the preceding are hereby excluded, without prejudice to any statutory or other contractual claims. This applies in particular to any claims for compensation for damage, consequential damage or costs relating to the defect. The transport of the instrument to us and return transport during this period shall be at the customer's risk.

To assert a claim under this undertaking, the instrument involved in the complaint must be sent to us, at the customer's expense and risk, with a detailed description of the defect complained of, specifying the purchase order number, delivery note number and invoice number. Our address is as follows:
Wankelstrasse 60, 50996 Cologne, Germany.

Claims under this undertaking expire 6 months from assertion of the claim, but not before the end of the claim period. This undertaking and any claims, rights and obligations arising therefrom are subject solely to the material law of the Federal Republic of Germany, to the exclusion of standards of International Private Law and the UN Law on the Sale of Goods. Our general supply and payment terms apply as supplementary conditions.

Any serious incident that has occurred in relation to the device must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

5.8 Declaration of Conformity

The current Declaration of Conformity can be found on our website: <http://www.pfmmedical.com/downloads>

Part 6 – Disposal of the pfm Rotary 3006 EM when put out of use

Return and disposal of waste electrical equipment

The WEEE directive 2002/96/EC (WEEE = Waste Electrical and Electronic Equipment) was issued in Europe with the aim of reducing the growing quantity of electronic scrap from discarded electrical and electronic devices and disposing of it in an environmentally friendly manner. In Germany, the directive was integrated in the 2005 Electrical and Electronic Equipment Act (ElektroG), which implemented both the WEEE directive and the RoHS (Restriction of Hazardous Substances in Electrical and Electronic Equipment).

Histotechnology-devices

As a registered manufacturer (WEEE Reg. No. DE85819911), we take back old histological electrical devices, manufactured by pfm medical gmbh, after they have been decommissioned and ensure that they are disposed of properly. For disposal, we work together with the certified specialist disposal company enretec GmbH, whose disposal facility is also registered with the ear foundation.

Arrange the return

To arrange for the return of your electrical device, please proceed as follows:
Contact us to arrange the return

- ▶ by phone: +49 2236 9641-0
- ▶ by e-mail: service[at]pfmmedical.com
- ▶ by mail: pfm medical gmbh, Abteilung Histotechnology, Wankelstraße 60, 50996 Köln, Germany

Prepare the device for transport

Please use only the original packaging and instruct a carrier of your choice to return the device to the following address:

pfm medical gmbh
Abteilung Histotechnology
Wankelstraße 60
50996 Köln
Germany

Costs

The costs of disposal are borne by us as the manufacturer. Unless otherwise agreed, the transport and packaging costs are to be covered by the owner/user of the device.

WARNING

Caution! If the device has come into contact with infectious material and cannot be disinfected or decontaminated, we cannot accept it for disposal.

Part 7 – Accessories



Universal cassette clamp

REF 010003

- ▶ Object clamp to be used with orientable adapter or fixed adapter
- ▶ Quick-release clamp system for all available cassettes
- ▶ Crossways and lengthways clamping possible



Universal cassette clamp, orientable

REF 013019

- ▶ Orientable version with adapter
- ▶ Quick-release clamp system for all available cassettes
- ▶ Crossways and lengthways clamping possible



Standard object clamp

REF 010004

- ▶ Object clamp to be used with orientable adapter or fixed adapter
- ▶ Maximum sample size: 40 x 45 mm



Standard object clamp, orientable

REF 010044

- ▶ Orientable version with adapter
- ▶ Maximum sample size: 40 x 45 mm



Universal cassette clamp large

REF 013022

- ▶ Quick-release clamp system for large embedding cassettes
- ▶ Crossways and lengthways clamping possible
- ▶ Object clamp to be used with orientable adapter or fixed adapter



Universal cassette clamp large, orientable

REF 013020

- ▶ Quick-release clamp system for large embedding cassettes
- ▶ Crossways and lengthways clamping possible
- ▶ Orientable version with adapter



Universal cassette clamp large slim

REF 013024

- ▶ Quick-release clamp system for super mega slim embedding cassettes
- ▶ Crossways and lengthways clamping possible
- ▶ Non-orientable version for direct adaptation



Universal cassette clamp large slim, orientable

REF 013023

- ▶ Quick-release clamp system for super mega slim embedding cassettes
- ▶ Crossways and lengthways clamping possible
- ▶ Orientable version with adapter



Standard object clamp large

REF 013010

- ▶ Non-orientable version for direct adaptation
- ▶ For larger samples
- ▶ Maximum sample size: 45 x 60 mm



Standard object clamp large, orientable

REF 013009

- ▶ Orientable version with adapter
- ▶ For larger samples
- ▶ Maximum sample size: 45 x 60 mm



Adapter orientable

REF 013016

- ▶ The orientable adapter is suitable for pfm Slide and Rotary microtomes
- ▶ For object orientation



Adapter fixed

REF 030055

- ▶ The fixed adapter is used for direct adaptation



Distance adapter object orientation 6 mm

REF 010046

- ▶ Suitable for all pfmmedical object clamps
- ▶ Distance 6 mm



Foil clamp, non-orientable

REF 030056

- ▶ Spring loaded jaws secure the sample in position
- ▶ To be used with orientable adapter or fixed adapter



Foil clamp, orientable

REF 030057

- ▶ Spring loaded jaws secure the sample in position
- ▶ Orientable version with adapter



Round sample clamp, orientable

REF 030058

- ▶ For sectioning of round samples (6 mm, 15 mm and 25 mm diameter)
- ▶ Orientable version with adapter
- ▶ Without inserts



Insert for round sample clamp, 6 mm

REF 030059

- ▶ Diameter 6 mm



Insert for round sample clamp, 15 mm

REF 030060

- ▶ Diameter 15 mm



Insert for round sample clamp, 25 mm

REF 030061

- ▶ Diameter 25 mm



V-insert

REF 030062

- ▶ For round samples of various diameters
- ▶ To be fitted in the standard object clamp



R-insert

REF 030063

- ▶ For round samples with a diameter of 24-25 mm
- ▶ To be fitted in the standard object clamp



Standard knife holder, upper jaw

REF 030064

- ▶ For use with steel and tungsten carbide knives up to 16 cm length
- ▶ With finger guard



Disposable blade holder, upper jaw, low profile

REF 030052

- ▶ For use with low profile microtome blades
- ▶ Removable pressure plate
- ▶ Easy to clean



Disposable blade holder, upper jaw, low profile black

REF 030072

- ▶ For use with low profile microtome blades
- ▶ Removable black pressure plate
- ▶ Easy to clean



Disposable blade holder, upper jaw, high profile

REF 030053

- ▶ For use with high profile microtome blades
- ▶ Removable pressure plate
- ▶ Easy to clean



Disposable blade holder, upper jaw HP TC

REF 030065

- ▶ For the use of Feather® high profile microtome blades (wide blades) with tungsten carbide blade edge
- ▶ Removable pressure plate
- ▶ Suitable for the pfm Rotary 3006 EM
- ▶ Easy to clean



Disposable blade holder, upper jaw, low profile lateral

REF 030071

- ▶ The upper jaw of the disposable blade holder can slide sideways (laterally)
- ▶ Three defined stop positions (left, middle, right)
- ▶ For use with low profile microtome blades
- ▶ Removable pressure plate
- ▶ Easy to clean



Pressure plate black

REF 030069

- ▶ Material: Black anodised aluminium
- ▶ Suitable for disposable blade holder 030052 and 030071



Blade holder base

REF 030051

- ▶ Designed to hold the disposable blade holder upper jaw
- ▶ Fixed to the microtome base plate using the knife holder mounting device
- ▶ Removable
- ▶ Easy to clean



Field magnifier set with illumination

REF 030066

- ▶ Magnification 2.5 x
- ▶ To facilitate the knife/microtome blade alignment and approach before sectioning
- ▶ Optimal view position



Foot switch

REF 030054

- ▶ Suitable for the pfm Rotary 3006 EM



Section waste tray

REF 030050

- ▶ Material: Anodised aluminium
- ▶ Easy to clean



Control panel

REF 030068

- ▶ Control panel for the pfm Rotary 3006 EM

For your notes

Contact

Do you have any questions?
Our Customer Solutions Team will be happy to advise you.

 service@pfmmedical.com

 +49 (0)2236 9641-220

 +49 (0)2236 9641-51



pfm medical gmbh
Wankelstraße 60
50996 Köln
Germany

