Implantable vascular access port

en	Instructions for Use	5
de	Gebrauchsanweisung	17
fr	Notice d'utilisation	29
es	Instrucciones de uso	41
it	Istruzioni per l'uso	53
nl	Gebruiksaanwijzing	65
sv	Bruksanvisning	77
el	Οδηγίες χρήσης	89
tr	Kullanma Talimatları	101
cs	Návod k použití	113
hu	Használati útmutató	125
ru	Инструкция по применению	137
қаз	Пайдалану туралы нұсқаулық	150
pl	Instrukcja użytkowania	163

Explanation of the symbols on label and packaging // Erklärung der Symbole auf Etikett und Verpackung // Explication des symboles sur l'étiquette et l'emballage // Explicación de los símbolos en la etiqueta y el envoltorio // Illustrazione dei simboli riportati su etichetta e confezione // Verklaring van de symbolen op het etiket en de verpakking // Förklaring av symboler på etiketten och förpackningen // Εξήγηση των συμβόλων πάνω στην ετικέτα και τη συσκευασία // Etiket ve ambalaj üzerindeki sembollerin açıklamaları // Vysvětlivky k symbolům na etiketě a balení // A címkén és a csomagoláson lévő szimbólumok magyarázata // Объяснение символов на этикетке и упаковке // Затбелгі мен бумадағы таңбалардың түсініктемелері // Objaśnienie symboli na etykietcie i opakowaniu

Reference number // Artikelnummer // Référence catalogue // Número catalogue // Numero di catalogo // Artikelnummer // Bokningsnummer // Αριθμός παραγγελίας // Sipariş numarası // Objednací číslo // Rendelési szám // Номер заказа // Тапсырыс нөмірі // Numer katalogowy

Lot number // Charge // Code de lot // Código de lote // Codice del lotto // Chargennummer // Parti- nummer // Αριθμός παρτίδας // Parti numarası // Číslo šarže // Gyártási tétel száma // Номер партии // Топтама нөмірі // Kod partii

Read Instructions for Use carefully // Gebrauchsanweisung beachten // Consulter les instructions d'utilisation // Consultense las instrucciones de uso // Consultare le istruzioni per l'uso // Volg de instructies in de gebruiksaanwijzing // Läs bruksanvisningarna noggrant // Διαβάστε προσεκτικά τις οδηγίες χρήσης // Kullanma kılavuzunu dikkate alınız // Respektujte Návod k použití // Vegye figyelembe a használati útmutatót // Соблюдать инструкцию по применению // Пайдалану бойынша нұсқаулықты сақтау // Zájrzij do instrukcji używania

Protect from direct sunlight // Vor Sonnenlicht schützen // Conserver à l'abri de la lumière du soleil // Manténgase fuera de la luz del sol // Tenere lontano dalla luce // Product niet blootstellen aan direct zonlicht // Skydda från direkt solljus // Να μην εκτεθούν στο άμεσο ηλιακό φως // Doğrudan güneş ışığından koruyunuz // Chraňte před přímým slunečním zářením // Óvja a közvetlen napfénytől // Не допускать попадания прямых солнечных лучей // Тікелей күн сәулесінен қорғау // Trzymać z dala od światła słonecznego

Store in a dry place // Trocken aufbewahren // Craint l'humidité // Manténgase seco / Mantenere asciutto // Droog bewaren // Förvara torr // Να διατηρείται σε ξηρό μέρος // Kuru ortamda saklayınız // Uchovávejte v suchu // Száraz helyen tartandó // Хранить в сухом месте // Құрғақ жерде сақтау // Chronić przed wilgocią

Expiry Date // Verwendbar bis // Date limite d'utilisation // Fecha de caducidad // Data di scadenza // Vervaldatum // Sista förbrukningsdag // Ημερομηνία λήξεως // Son kullanma tarihi // Použitelné do // Lejárati idő // Срок хранения // Сақтау мерзімі // Użyć do daty

For single use only // Nicht wiederverwenden // Ne pas réutiliser // No reutilizar // Solo per impiego monouso // Alleen bestemd voor eenmalig gebruik // Använd endast en gång // Για μία χρήση μόνο // Bir defalık kullanım içindir // Určeno k jednorázovému použití // Kizárólag egyszeri használatra alkalmas // Только для одноразового использования // Бір жолғы пайдалануға ғана // Nie używać powtórnie

Sterilized by Ethylene Oxide // Sterilisiert mit Ethylenoxid // Stérilisé avec de l'oxyde d'éthylène // Esterilizado utilizando óxido de etileno // Sterilizzato con ossido di etilene // Gesteriliseerd met ethyleenoxide // Steriliserat med etylenoxid // Αποστειρωμένο με αιθυλενοξείδιο // Etilen oksit ile sterilize edilmiştir // Sterilizováno ethylenoxidem // Etilénoxiddal sterilizált // Простерилизовано оксидом этилена // Этилен оксидімен зарарсыздандырылған // Sterylizowane tlenkiem etylenu

Do not resterilise // Nicht erneut sterilisieren // Ne pas restériliser // No esterilizar // Non risterilizzare // Niet opnieuw steriliseren // Omsteriliserá inte // Μην αποστειρώνετε εκ νέου // Tekrar sterilize etmeyin // Nepoužívejte nesterilizujte // Ne sterilizálja // Не подвергать рестерилизации // Қайта зарарсыздандыруға болмайды // Nie resterylizować

Manufacturer // Hersteller // Fabricant // Fabricante // Fabbicante // Fabrikant // Tillverkare // Κατασκευαστής // Üretici // Výrobce // Gyártó // Изготовитель // Өндіруші // Wytwórca

Latex free // Latex-frei // Sans latex // Sin látex // Non contiene lattice // Latex-vrij // Latexfri // Χωρίς λάτεξ // Latexs içermez // Latex zdarma // Latex mentes // Не содержит латекса // Кұрамында латекс жоқ // Nie zawiera lateksu

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Implantable vascular access port

Product description

This instruction for use is valid for the following implantable vascular access ports:

T-Port model range	Jet Port® Plus model range	TitaJet™ light model range
T-Port	Jet Port® Plus	TitaJet™ light
T-Port Low Profile	Jet Port® Plus Low Profile	TitaJet™ light Low Profile Contrast
T-Port Contrast	Jet Port® Plus II Contrast	TitaJet™ light II Contrast
T-Port Low Profile Contrast		

Implantable vascular access ports are fully implantable, subcutaneously located and facilitate repeated vascular access intended for the delivery of medications, intravenous fluids, parenteral nutrition solutions or blood products and blood sampling through the use of a Huber needle (specific non-coring needle). The implantable vascular access port can be implemented for outpatient treatment. It remains in place throughout the duration of patient treatment and can be explanted when needed.

The implantable vascular access port consists of a port chamber ①, a septum ② and an outlet tube ③ where a catheter ④ is secured by a radiopaque click connector ⑤. Refer to Figure 1: Illustration of the implantable vascular access port.


Figure 1: Illustration of the implantable vascular access port.

The materials of the different components are as follows:

Components					
Model range	Port body ①	Septum ②	Outlet tube ③	Catheter ④	Radiopaque click connector ⑤
T-Port	Titanium	Silicone	Stainless steel	Polyurethane or silicone	Radiopaque polyurethane resin
Jet Port® Plus	Polyoxymethylene				
TitaJet™ light	POM with a titanium base plate (insert) inside the port chamber				


Table 1: Component materials

The patient population is not limited to a specific age, sex or weight. Only the adequacy of vascular anatomy is pivotal and needs to be evaluated by the physician individually.

Drilled holes are provided at the edge of the port chamber's base plate to secure the port chamber in subcutaneous tissue by sutures. The septum has a high pressure stability and holds the needle that has punctured the septum securely in position. A high number of punctures is possible using a Huber needle (e.g. JetCan™). The connectable catheters can be shortened to the required length at the proximal end (closest to the port). Various versions, sizes and materials are available for the catheters to address different medical issues. The silicone and polyurethane catheters, as well as the click connector, are fully radiopaque.

Note: Only Huber needles (e.g. JetCan™) should be used to puncture the septum. The special design of the needle tip prevents silicone particles from being punched out.

The implantable vascular access ports include the following components:

1x Patient ID card

1x Set of serial number label

1x Instruction for use

The implantable vascular access contrast ports, which are power-injection compatible, have a purple port body outer surface and a CT radiopaque indicator at their back. A purple patient bracelet is provided for an easy identification of the contrast port after implantation.

General Information and Warnings

These instructions for use and the information on the packaging should be read carefully before each use.

- The instructions for use is intended for the attending physician and medical staff. The patient ID card booklet is intended for patients, their family members, hospital staff and also serves as general information.
- The patient ID card must be filled out by the attending physician. The patient should always carry the patient ID card to ensure that all the necessary information is quickly available in the event of complications.
- A physician with appropriate knowhow is responsible for evaluating the physical condition and suitability of the patient as well as the condition of the implantable vascular access port and its position in the vessel in relation to the administration of contrast medium and flow rate.
- All materials used are biocompatible.

WARNINGS

- ☞ **Do not use the product if the packaging has been opened or is damaged, if you are not sure that it is sterile, or if the expiry date has passed.**
- ☞ **Each product is packed separately and is delivered in an EO-sterilised and pyrogen-free condition. It is intended for single use only. Do not reuse or implant a second time, reprocess or resterilize. Re-use, reprocessing or resterilization of single-use devices may result in degraded performance or a loss of functionality. Reuse of single-use devices may result in exposure to pathogens such as viruses, bacteria, fungi or prions.**

Sterility

It is essential to ensure sterile, strictly aseptic conditions for all procedures.

Avoidance of damage to system components and/or avoidance of injury to the patient during the implantation

- Ensure that there is no risk of mechanical damage to system components, particularly from sharp objects. Only rounded, atraumatic clamps and forceps should be used.
- Avoid perforating, tearing or fracturing the catheter when using the guidewire.
- Never use a catheter if there is any suspicion of mechanical damage.
- During the implantation, ensure that there is no sharp angle or kinking in the catheter when using the introducer system.
- Insert the introducer and the catheter in order to avoid penetration into the pleura of the lung.
- Do not allow the percutaneous introducer to remain in the blood vessel without the internal support of a catheter or dilator.
- Simultaneously advance the introducer and dilator with a rotational motion.

Indication

The fully implantable vascular access ports, subcutaneously located, are indicated for patient therapies requiring repeated vascular access.

The following substances can be administered:

- Medications (e.g. cytostatics)
- Infusion solutions (e.g. NaCl)
- Blood products
- Substances for parenteral nutrition
- Administration of contrast medium

In certain circumstances implantable vascular access ports are also suitable for:

- Taking blood samples

WARNINGS

 The indication of the administration of contrast medium (contrast ports, power injection compatible, with a catheter- dependent flow rate of up to 5 ml/s, but always under a maximum pressure of 300 psi (20.7 bar)) shows that the implantable vascular access ports can withstand this application (always refer to paragraph Technical Specifications, for suitability verification). This does not imply that the application is suitable for a particular patient, or for a particular infusion set.

Contraindications

Implantable vascular access ports may not be implanted in the case of the following contraindications:

- Inadequate physical condition that prevents optimum functioning of the system
- Bacteremia
- Sepsis
- Known or suspected allergy to the materials used
- Severe chronic obstructive pulmonary disease
- Prior or scheduled radiation therapy in the region of the implantation site
- Previous venous thrombosis
- Planned or already performed vascular surgical or interventional radiological procedure in the region of the implantation
- Local tissue factors that might prevent proper functioning of the system or operational access
- Severe coagulation disorders
- Patient's body size insufficient to accommodate the size of the implanted device

Possible Complications

Implantable vascular access ports facilitate vital access to the vascular system of critically ill patients. However, they can cause serious complications such as:

- | | |
|---|---|
| <ul style="list-style-type: none"> • Air embolism • Catheter embolism • Hemorrhaging • Haematoma formation • Brachial plexus injury • Cardiac arrhythmia • Catheter or port erosion through the skin • Catheter or port occlusion (e.g. Blood reflux dysfunction) • Catheter or port sepsis • Endocarditis • Rotation or extrusion of the system • Extravasation • Fibrin thrombi formation • Phlebitis (e.g. Thrombophlebitis) • Pneumothorax | <ul style="list-style-type: none"> • Hemothorax • Hydrothorax • Intolerance reaction to the implanted biomaterial • Inflammation, necrosis or scarring of skin over implant area • Lesion/perforation of vessels and soft tissue • Thoracic duct injury • Pericardial or cardiac tamponade (thrombus formation) • Spontaneous catheter tip malposition or retraction • Thromboembolism • Vascular thrombosis • Vessel erosion • Risks generally associated with local and general anesthesia, surgery and post-operative recovery |
|---|---|

Follow-up

After the implantation, implantable vascular access ports must be checked regularly upon their function and position. They should be used only by physicians and other medical staff who are familiar with the procedure. Medical staff should have sufficient experience of the possible complications involved. Complications can arise at any time during or after each treatment cycle.

Implantation Instructions

Sterile conditions must be strictly observed during the entire implantation procedure.

There are many different techniques for implanting vascular access ports and inserting catheters into the vascular system. Techniques include for example the percutaneous procedure (Seldinger technique), the introducer technique or the catheter tunneling procedure. In the following the Surgical (open) implantation procedure will be described. The choice of technique is the responsibility of the physician performing the implantation.

Note: A thrombosis prophylaxis is required for arterial vessel access, e.g. via the subclavian artery or femoral artery.

WARNINGS

 **Review of literature shows that insertion of implantable vascular access ports on peripheral vein increases the risk of venous sclerosis.**

Surgical (open) implantation procedure

- ❶ Document the implantation procedure in full, including the production number, re-order number and batch number.
- ❷ Decide on the location and side of body for the port implantation (standard method: right infra-clavicular fossa). However, a different site may be more suitable because of individual patient factors. The selected port site in a specific anatomical area or part of the body should ensure port stability that does not interfere with patient mobility (such as the wearing of particular items of clothing, e.g. a bra). An additional factor to remember is the thickness of the skin or subcutaneous adipose tissue over the port chamber. A tissue thickness of 0.5 - 2.0 cm should be preferred. If the skin or subcutaneous tissue is too thin, it can lead to erosion above the port. If the skin or subcutaneous tissue is too thick, access to the port chamber can become extremely difficult.
- ❸ Prepare a sterile operating field.
- ❹ Perform appropriate anesthesia.
- ❺ Prepare the catheter and remove the air from the catheter, connect the rinsing needle to the catheter at the end nearest the port. Rinse the catheter with heparinized saline solution and clamp the catheter at the port end.
Note: The catheter should be clamped on a segment that will be cut off and discarded before the catheter is connected to the port chamber.
- ❻ Prepare and expose the selected vessel. Following exteriorization and preparation, ensure there is no bleeding and/or air aspiration. Only then should an incision be made in the vessel, either directly (terminally), or with a small pocket incision. If a catheter is to be introduced into the subclavian vein, refer to the warnings on catheter compression in the pinch-off area.
- ❼ Introduce the appropriate end of the vein lifter through the incision into the vessel.
- ❽ Advance the catheter along the underside of the vein lifter into the vessel.
- ❾ Prepare the subcutaneous pocket for the port implantation. Part of the pocket can be prepared bluntly. Ensure that the pocket is large enough to receive the port. The port should not be positioned too close to the vessel incision.
- ❿ Advance the catheter into the vessel as far as the selected infusion site. The tip of the catheter should come in the vena cava superior. A suitable technique should be employed to verify that the catheter tip is in the correct position.
Note: If ligatures or other methods are used to secure the catheter on the vessel wall, care must be taken to avoid occlusion or damage to the catheter.
- ⓫ Withdraw the vein lifter.

- 12 Create a tunnel through the subcutaneous tissue from the port pocket to the site of the vessel incision. Withdraw the catheter end in the direction of the port pocket. Cut the catheter to the desired length at an approx. 90° angle to the catheter lumen. This is the only way to ensure the port can be securely connected at a later stage. Place the click connector locking mechanism over the catheter.
- 13 Before connecting, remove the air from the port chamber and rinse the port chamber according to the facility protocols (for instance: heparinized saline solution) using a Huber needle to puncture the septum.
- 14 Rinse all system components with a suitable rinsing solution.
- 15 Connect the catheter to the port (refer to paragraph Connecting the Catheter to the Port).
- 16 Place the port in the subcutaneous pocket at a suitable distance from the incision site and secure it to the underlying fascia with one non-resorbable, monofilament suture per suture hole. Ensure that the catheter is implanted loosely enough to allow slight movement and that there are no kinks in it.
- 17 After suturing the port in the pocket, flush the wound with a suitable antibiotic solution.
- 18 Check the flow in the catheter using a Huber needle and a syringe to ensure that the flow is not obstructed, that there is no leakage and that the catheter is correctly positioned. Aspirate to confirm that blood can be withdrawn.
- 19 Flush and heparin lock the implantable vascular access port.
- 20 Close the incision site.

Note: If the catheter and the click connector are connected and then disconnected, the port end of the catheter should again be shortened at an approx. 90° angle to the catheter lumen to ensure a secure connection.

WARNINGS

- ☞ **It is important to ensure that the catheter is correctly positioned before attaching the catheter lock. An incorrectly positioned catheter can lead to a system malfunction (e.g. catheter disconnection or extravasation).**
- ☞ **The catheter should be straight or slightly curved and sit at the vessel incision without any pressure or tension (even when the patient moves). If there are kinks, the catheter should be straightened and if necessary the kinked segment should be removed. The catheter may be damaged if the catheter lock is placed over a kinked catheter.**

Connecting the Catheter to the Port

- 1 Place the radiopaque click connector locking mechanism over the catheter.
- 2 Before connecting, remove the air from the port chamber and rinse the port chamber with heparinized saline solution using a Huber needle to puncture the septum.
- 3 Rinse all system components with a suitable rinsing solution.
- 4 Pull click over catheter as shown for orientation.
- 5 Advance the catheter in a straight line over the pin attached to the port housing.
- 6 Advance the catheter lock that was previously placed over the catheter until it clicks into place.


Use Instructions

- ❶ Sterile conditions must be maintained during the entire procedure.
- ❷ Check that the catheter and port chamber are securely connected.
- ❸ Use only Huber needles, such as JetCan™, to puncture the port septum.

Note: The Huber needle must not be left open in the system (air aspiration).

- ❹ JetCan™ Huber needles are designed and authorized for single use only.
- ❺ No syringes smaller than 10 ccm (ml) should be used. If occlusion of the system occurs, excessive pressure can damage the port septum or the catheter.
- ❻ Before each treatment cycle, check that the port chamber is in the correct position and for any signs of wound infection or extravasation in the vicinity of the puncture site (refer to paragraph Indications for Explanation).

- ❼ For the infusion of medication, refer to the individual pharmaceutical manufacturer's instructions.
- ❽ All punctures of the implantable vascular access ports must be carried out under sterile conditions.
- ❾ After palpating the outer ring of the port chamber, a Huber needle, such as the JetCan™, should be used to puncture the port chamber at a slight angle through the septum.
- ❿ If more than one drug is used during a therapy cycle, the system should be rinsed between applications with 5 - 10 ccm (ml) physiological saline solution to avoid any interaction between the drugs.
- ⓫ After each infusion, injection or bolus application, the system must be rinsed with 5 - 10 ccm (ml) of physiological saline solution and then locked with a solution according to the facility protocols: for instance a heparin lock with heparinized saline solution. This will prevent occlusion of the catheter.
- ⓬ The Huber needle should be withdrawn under slight, positive pressure to reduce the risk of blood aspiration into the catheter lumen.

Application of fluids

The following instructions must be followed before any application of fluids:

- ❶ Only Huber needles, such as JetCan™, may be used for any access to the implantable vascular access port. Check that the implantable vascular access port is clear by flushing it with approx. 5 ccm (ml) NaCl 0.9 %. The flushing process must be smooth flowing. Next, a sample aspiration should be performed on the implantable vascular access port until venous blood is visible. After that, the implantable vascular access port should be rinsed with approx. 10 ccm (ml) NaCl 0.9 %. Normally, small thrombi are deposited at the catheter tip, which can act as a kind of valve mechanism. If it is not possible to aspirate blood, in the majority of cases it is still possible to apply a fluid through the implantable vascular access port without too much difficulty.
- ❷ If it is not possible to flush the system, the position of the Huber needle (JetCan™), the catheter or the entire system should be checked (see also catheter occlusion).
- ❸ A n application can be performed as a bolus injection or an infusion after the system has been successfully rinsed.
- ❹ A fter the application, the system must be rinsed and locked (refer to paragraph Maintenance Procedure).

It must be possible to inject fluid without resistance, since otherwise there might be a blockage. Under no circumstances should an attempt be made to rinse clear the occlusion under pressure, as this could result in damage to the port or the catheter. Should the tissue become swollen while fluids are being administered, terminate the procedure immediately (extravasation) and take appropriate steps.

Recommendation

The indwelling time of a Huber needle is up to 72 hours and when blood products and lipid emulsions are infused, the Huber needle must be changed after 24 hours. The indwelling time of a Huber needle could be defined according care center's applicable policies and procedures, under their own responsibility. Those policies and procedures should be based on patient's medical condition, laboratory tests, prior experience or other recognized medical evidences.

WARNINGS

-  **When the system is not in use, it must be completely filled (locked) at regular intervals according to the facility protocol: for instance a heparin lock (refer to paragraph Maintenance Procedure).**
-  **Verify needle length is correct based on port reservoir depth, tissue thickness, and the thickness of any dressing beneath the bend of the needle; if too long, needle and/or port may be damaged during insertion; if too short, needle may not completely pierce port septum, and medication may be delivered into surrounding tissue and/or needle may be blocked.**

Additional User Instructions: Application of Contrast Medium for Computer Tomography

- The port must be rinsed before the injection of contrast medium and it must operate properly (see above).
- Before the contrast medium is injected, we recommend checking by means of a topogram that the radiopaque port catheter is in the correct position in the vessel, assuming this cannot be determined from a recent chest X-ray.
- The contrast medium should be preheated to body temperature before being injected.
Note: If this is not done a malfunctioning of the implantable vascular access port could result.
- To ensure an adequate flow rate, we recommend the use of a 19 G (1.1 mm) Huber needle of the JetCan™ range.
- The settings at the infusion pump should be such that, depending on the viscosity of the contrast medium, a maximum pressure of 300 psi (20.7 bar) is not exceeded.
- Advise the patient to inform you immediately upon any sign of pain, or if he/she begins to feel unwell during the injection. If this happens, discontinue the injection at once.
- Following the use of contrast media, the system must be rinsed immediately with at least 20 ml NaCl 0.9 %.
- After the application, the system must be rinsed and locked (refer to paragraph Maintenance Procedure).

Additional User Instructions: Blood Transfusions or Aspiration

If the implantable vascular access port is used for blood transfusions or taking venous blood samples, the instructions for use basically apply. Additional, specific steps are the following:

- Aspirate 3 ccm (ml) blood and then discard the sample.
- Aspirate the quantity of blood required for the examination.
Note: Aspiration with Huber needles smaller than 22 G can result in inaccurate laboratory results.
- Immediately after taking the blood sample, the system should be thoroughly rinsed with 30 ccm (ml) heparinized saline solution. During the rinsing procedure, the injection needle should be rotated inside the port chamber.
- After the application, the system must be rinsed and locked (refer to paragraph Maintenance Procedure).

Note: Using the implantable vascular access port to take frequent blood samples considerably reduces the service life and functionality of the system.

Maintenance Procedure

To prevent clot formation and catheter blockage, implanted ports should be filled with sterile heparinized saline solution after each use. Please refer to paragraph Technical Specification: Table 2 for port priming volumes and Table 3 for catheter priming volumes to determine the volume required to lock the implantable vascular access port. When not using the implantable vascular access port, the heparin lock should be changed at least once every four weeks.

It is recommended in any event to rinse the implantable vascular access port with a suitable rinsing solution as described to ensure that the system remains clear.

Heparin lock

Recommended heparin concentration to produce a heparin lock for every implantable vascular access port after the rinsing procedure: 100 IU/ccm heparinized NaCl.

WARNINGS

 **Heparin must not be administered as a locking or rinsing solution if the patient suffers from heparin hypersensitivity or heparin-induced thrombocytopenia (HIT).**

Note: The maintenance of the implantable vascular access port should be executed according these maintenance recommendations or according care center's applicable policies and procedures, under their own responsibility. Those policies and procedures should be based on patient's medical condition, laboratory tests, prior experience or other recognized medical evidences.

Special Procedure in the Event of a Genuine Catheter Occlusion

In many instances, functional occlusions (compressions) of the catheter can be remedied by repositioning or moving the patient.

Micro-thrombus formation and drug residual are normally flushed out of the system during the rinse procedure. This is why no syringes with a capacity of less than 10 ccm (ml) should be used.

If serious catheter complications do arise and other sources can be excluded, the following procedures may help:

Fibrinolytic solution needs only be injected if the system is clotted with blood.

- Prepare a solution of urokinase and physiological saline solution with a concentration of 5000 U/ccm.
- A small amount of the solution should be injected using a 10 ccm (ml) syringe and a Huber needle. Rinse carefully, alternately aspirating and injecting solution.
- Wait 5 - 10 minutes and then repeat the process. It may already be possible either to aspirate or flush.
- If necessary, the procedure can be repeated 5 minutes later and at subsequent 5-minute intervals, until the catheter is clear again.
- After the application the system must be rinsed and locked (refer to paragraph Maintenance Procedure).

Note: The concentrations specified above are intended merely as a recommendation.

Depending on the coagulation situation of the patient, it is the attending physician's responsibility to adjust the dosage and administration of the fibrinolytic solution used, referring to the fibrinolytic manufacturers' instructions. Never attempt to rinse the system clear with force. This can lead to serious complications.

Indications for Explantation of the Implantable Vascular Access Port

Observance of the above instructions should ensure the long life of the implantable vascular access port. After five years, the elastic properties of the septum may diminish.

Below is a non-exhaustive list of indications that may necessitate explantation of the implantable vascular access port:

- Fever due to infection (e.g. bacteremia) with no evidence of another focus
- Haematoma, Inflammation, necrosis or scarring of skin over implant area
- Existing endoallergic reaction when using the implantable vascular access port
- Suspected endocarditis
- Suspected embolism
- Definite, non-reversible occlusion of the implantable vascular access port

Duration of Implantation of the Vascular Access Port

The port and catheter system should be removed at the end of treatment.

The duration of the vascular access port is mostly dependant on the catheter. Fracture is the most common complication of ageing. Risk of fracture may increase over time depending on the material of the catheter, the pathway and the entry point of the catheter into the vessel.

Removal of the Vascular Access Port

When removing the vascular access port, care must be taken not to fracture the catheter. If the catheter is sutured to the vessel, the sutures should first be removed. Control the catheter while removing the vascular access port from the port pocket. Be vigilant if there is excessive resistance to the removal of the catheter. Catheters may become encapsulated and attached to the vein wall. Should this occur and it is not possible to remove the catheter without risking catheter fracture, or if the catheter is fractured, the advice of an Interventional Radiologist, Surgeon or other Physician with endo-luminal experience should be sought.

Technical Specifications

Please refer to Table 3 for port priming volumes and Table 4 for catheter priming volumes to determine the volume required to lock the implantable vascular access port.

Model range	Port priming volume in ml (PORT+PIN)
Jet Port® Plus Low Profile	0.4
Jet Port® Plus	0.6
Jet Port® Plus II Contrast	0.8
TitaJet™ light Low Profile Contrast	0.4
TitaJet™ light	0.5
TitaJet™ light II Contrast	0.7
T-Port Low Profile (Contrast)	0.4
T-Port (Contrast)	0.8

Table 2: Port priming volume

Catheter material	Catheter size in French	Internal diameter in mm	Outer diameter in mm	Priming volume in ml/10 cm
Polyurethane	4.8	1.0	1.6	0.1
Polyurethane	6.6	1.3	2.2	0.15
Silicone	6.6	1.0	2.2	0.1
Polyurethane	8.0	1.6	2.6	0.2
Silicone	8.4	1.6	2.8	0.2

Table 3: Catheter priming volume

Table 4 displays the gravity flow rates of the different catheters for the implantable vascular access ports.

Catheter material	Catheter size in French	Gravity flow rate in ml/min
Polyurethane	4.8	≤ 13*
Polyurethane	6.6	≤ 13*
Silicone	6.6	≤ 22*
Polyurethane	8.0	≤ 22*
Silicone	8.4	≤ 25*

Table 4: Catheter gravity flow rate

*Flow rate for 1 m of water column ($\gamma = 98.1$ mbar), Huber needle 19 G x 25 mm - catheter length 40 cm - according to ISO 10555-6:2015.

Table 5 displays the power injection flow rates of the different implantable vascular access ports.

Model range	Contrast/ non-Contrast	Power injection flow rate in ml/s
Jet Port® Plus II Contrast	Contrast	≤ 5 ml/s**
TitaJet™ light II Contrast	Contrast	≤ 5 ml/s**
TitaJet™ light Low Profile Contrast	Contrast	≤ 5 ml/s**
TitaJet™ light Low Profile F 4.8	Contrast	≤ 3 ml/s**
T-Port (Contrast)	Contrast	≤ 5 ml/s**
T-Port Low Profile (Contrast)	Contrast	≤ 5 ml/s**

Table 5: Power injection flow rate

**Max flow rate at a pressure under 300 psi (20.7 bar) for a liquid injected at a viscosity of 6.1 Cp at 37 °C, Huber needle 19 G x 20 mm - catheter length 40 cm - according to ISO 10555-6:2015.

Table 6 displays the recommended maximum number of punctures for the different implantable vascular access ports.

Model range	Recommended maximum punctures
Jet Port® Plus Low Profile	≤ 120***
Jet Port® Plus	≤ 780***
Jet Port® Plus II Contrast	≤ 780***
TitaJet™ light Low Profile Contrast	≤ 120***
TitaJet™ light	≤ 780***
TitaJet™ light II Contrast	≤ 780***
T-Port Low Profile (Contrast)	≤ 200***
T-Port (Contrast)	≤ 780***

Table 6: Maximum number of punctures recommended per model range

*** The model range is compliant to ISO Standard 10555-6 §4.5.3, Annex D, no leakage after for 1000 punctures per cm² of septum surface area, Huber needle 19 G.

MRI Compatibility Information

Through non-clinical testing with a 3 Tesla (T) MR System it has been shown that the pfm medical cpp implantable vascular access port is MR conditional (MR safe under specific conditions). A patient with this implant can be scanned safely immediately after placement under the following conditions:

- Static Magnetic field of 3 T or less
- Spatial gradient magnetic field of 720 Gauss/cm or less
- Maximum MR system reported whole body average specific absorption rate (SAR) of 3 W/kg for 15 minutes of scanning.

In non-clinical testing, the pfm medical cpp implantable vascular access port produced a temperature rise of 0.7 °C at a maximum. MR system-reported whole body average specific absorption rate (SAR) of 3 W/kg for 15 minutes of MR scanning in a 3 T MR system using a transmit/receive body coil (Excite, Software G3.0052B, General Electric Healthcare, Milwaukee, WI).

MR image quality may be compromised if the area of interest is in the exact same area or relatively closed to the position of the port. Therefore, optimization of MR imaging parameters to compensate for the presence of this implant may be necessary.

WARNINGS

-  **The pfm medical cpp implantable vascular access port is not tested on MRI systems operating at more than 3 T.**

Disposal after Use

After use, medical products and accessories pose a potential biological hazard. For this reason, the products and their accessories should be handled and disposed of in accordance with recognised medical procedure and in compliance with the relevant legal regulations and local ordinances.

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pfm medical cpp S.A.