Titanized Mesh Implants TiLOOP® Bra

de	Gebrauchsanweisung	8
en	Instructions for use	13
fr	Notice d'utilisation	17
es	Instrucciones de uso	22
it	Istruzioni per l'uso	26
nl	Gebruiksaanwijzing	31
da	Brugsanvisning	36
el	Οδηγίες χρήσεως	41
pt	Instruções de utilização	46
fi	Käyttöohjeet	51
sv	Bruksanvisning	56
et	Kasutuskorraldus	61
lv	Lietošanas instrukcija	66
lt	Naudojimo instrukcija	71
pl	Instrukcja użytkowania	76
sk	Návod na použitie	81
sl	Navodila za uporabo	86
cs	Návod k použití	91
hu	Használati utasítás	96
bg	Инструкция за употреба	101
ro	Instrucțiuni de utilizare	106
hr	Uputa o načinu upotrebe	111
no	Bruksanvisning	116
sr	Uputstvo za upotrebu	121
ru	ИНСТРУКЦИЯ ПО ИСПОЛЬЗОВАНИЮ	125
tr	Kullanma kılavuzu	130

de Erklärung der Symbole auf Etikett und Verpackung **en** Explanation of the symbols on label and packaging **fr** Explication des symboles sur l'étiquette et l'emballage **es** Explicación de los símbolos en la etiqueta y el envoltorio **it** Illustrazione die simboli riportati su etichetta e confezione **nl** Verklaring van de symbolen op het etiket en de verpakking **da** Forklaring af symbolerne på etiket og emballage **el** Εξήγηση των συμβόλων πάνω στην ετικέτα και τη συσκευασία **pt** Explicação dos símbolos na etiqueta e embalagem **fi** Etiketin ja pakkauksen symbolien selitys **sv** Förklaring av symboler på etiketten och förpackningen **et** Etiketil ja pakendil olevate sümbolite selgitus **lv** Simbolu uz etiķetes un iepakojuma izskaidrojums **lt** Simbolių ant etiketės ir pakuotės paaiškinimas **pl** objaśnienie symboli na etykiecie i opakowaniu **sk** Vysvetlivky symbolov na štítkoch a obaloch **sl** Pojasnila simbolov na etiketi in embalaži **cs** Vysvětlivky k symbolům na etiketě a balení **hu** Jelmagyarázat a címkehez és csomagoláshoz **bg** Обяснение на символите върху етикета и опаковката **ro** Explicația simbolurilor de pe etichetă și de pe ambalaj **hr** Objašnjenje simbola na etiketi i pakiranju **no** Forklaring av symboler på etikett og pakning **sr** Objašnjenje simbola na etiketi i na pakovanju **ru** Объяснение символов на этикетке и упаковке **tr** Etiket ve ambalaj üzerindeki sembollerin açıklamaları

de Medizinprodukt **en** Medical device **fr** Dispositif médical **es** Producto médico **it** Dispositivo medico **nl** Medisch product **da** Medicinsk produkt **el** Ιατρικό προϊόν **pt** Dispositivo médico **fi** Lääkinnällinen laite **sv** Medicinsk produkt **et** Meditsiinitoode **lv** Medicīniskā ierīce **lt** Medicinos produktas **pl** Produkt medyczny **sk** Zdravotnícky produkt **sl** Medicinski proizvod **cs** Zdravotnický výrobek **hu** Orvostechnikai eszköz **bg** Медицински продукт **ro** Dispozitiv medical **hr** Medicinski proizvod **no** Medisinsk produkt **sr** Medicinski proizvod **ru** Медицинский продукт **tr** Tibbi ürün

de Hersteller **en** Manufacturer **fr** Fabricant **es** Fabricante **it** Produttore **nl** Fabrikant **da** Producent **el** κατασκευαστής **pt** Fabricante **fi** Valmistaja **sv** Tillverkare **et** Tootja **lv** Izgatavotājs **lt** Gamintojas **pl** Producent **sk** Výrobce **sl** Proizvajalec **cs** Výrobce **hu** Gyártó **bg** Производител **ro** Producător **hr** Proizvođač **no** Produsent **sr** Proizvođač **ru** Изготовитель **tr** Üretici

de Herstellungsdatum **en** Date of manufacture **fr** Date de fabrication **es** Fecha de fabricación **it** Data di produzione **nl** Productiedatum **da** Fabrikationsdato **el** Ημερομηνία παραγωγής **pt** Data de fabrico **fi** Valmistuspäivä **sv** Tillverkningsdatum **et** Valmistamiskuupäev **lv** Izgatavošanas datums **lt** Pagaminimo data **pl** Data produkcji **sk** Dátum výroby **sl** Datum proizvodnje **cs** Datum výroby **hu** Gyártási idő **bg** Дата на производство **ro** Data fabricației **hr** Datum proizvodnje **no** Produksjonsdato **sr** Datum proizvodnje **ru** Дата изготовления **tr** Üretim tarihi

de Bestellnummer **en** Reference number **fr** N° de commande **es** Número de referencia para pedidos **it** Numero d'ordine **nl** Bestelnummer **da** Bestillingsnummer **el** Αριθμός καταλόγου **pt** Número de encomenda **fi** Viitenumero **sv** Beställningsnummer **et** Tellimisnumber **lv** Pasūtījuma numurs **lt** Užsakymo numeris **pl** Numer zamówienia **sk** Objednávacie číslo **sl** Številka naročila **cs** Objednací číslo **hu** Rendelési szám **bg** Номер на поръчка **ro** Număr de comandă **hr** Broj narudžbe **no** Bestillingsnummer **sr** Broj narudbe **ru** Номер для заказа **tr** Sipariş numarası

LOT

de Chargen-Nummer **en** Batch number **fr** N° de lot **es** Número de lote **it** Numero del lotto **nl** Chargennummer **da** Lotnummer **el** Κωδικός παρτίδας **pt** Número do lote **fi** Eränumero **sv** Satsnummer **et** Partii number **lv** Partijas numurs **lt** Partijos numeris **pl** Kod partii **sk** Číslo šarže **sl** Številka šarže **cs** Číslo šarže **hu** Gyártási tételszám **bg** Партиден номер **ro** Cod de lot **hr** Broj šarže **no** Produktpartinummer **sr** Serijski broj **ru** Номер партии **tr** Yük numarası

de Verwendbar bis **en** Use by **fr** Date limite d'utilisation **es** Úsese hasta **it** Da usare fino a **nl** Niet te gebruiken na **da** Anvendes senest **el** Ημερομηνία λήξης **pt** Utilizar de preferência antes de **fi** Viimeinen käyttöpäivämäärä **sv** Används före **et** Kasutatav kuni **lv** Izlietot līdz **lt** Galima naudoti iki **pl** Zużyć przed **sk** Použiteľné do **sl** Uporabno do **cs** Použitelné do **hu** Felhasználhatóság: **bg** Годно до **ro** A se utiliza până la **hr** Upotrebljivo do **no** Kan brukes inntil **sr** Upotrebljivo do **ru** Годен до **tr** Kullanma süresi

de Gebrauchsanweisung beachten **en** Follow Instructions for Use **fr** Respecter les instructions du mode d'emploi **es** Respete las instrucciones **it** Seguire le modalità d'uso **nl** Gebruiksaanwijzing raadplegen **da** Følg brugsanvisningen **el** Προσέξτε τις οδηγίες χρήσης **pt** Cumprir as instruções de utilização **fi** Katso käyttöohje **sv** Följ bruksanvisningen **et** Järgige kasutuskorraldust **lv** Ievērojiet lietošanas instrukcijas norādījumus **lt** Laikykites instrukcijų **pl** Zapoznaj się z instrukcją stosowania **sk** Dodržujte návod na použitie **sl** Upošteevajte navodila za uporabo **cs** Dodržujte návod k použití **hu** Vegye figyelembe a használati utasítást **bg** Обърнете внимание на инструкцията за употреба **ro** Respectați instrucțiunile de utilizare **hr** Poštivati uputu o uporabi **no** Overhold bruksanvisningen **sr** Obratiti pažnju na uputstva za upotrebu **ru** Соблюдайте инструкцию по эксплуатации **tr** Kullanım talimatına dikkat

de Nicht zur Wiederverwendung **en** Do not reuse **fr** Produit non réutilisable **es** No reutilizable **it** Non riciclabile **nl** Niet bestemd voor hergebruik **da** Må ikke genanvendes **el** Μόνο για μία χρήση **pt** Não reutilizável **fi** Ei saa käyttää uudelleen **sv** För engångsbruk **et** Mitte korduvkasutamiseks **lv** Nelietot atkārtoti **lt** Neskirta pakartotiniam naudojimui **pl** Nie używać powtórnie **sk** Nie k opätovnému použitiu **sl** Ni za ponovno uporabo **cs** Ne k opakovanému použití **hu** Nem újra felhasználható **bg** Не e за многократна употреба **ro** A nu se reutiliza **hr** Neprikladno za ponovnu uporabu **no** Ikke bruk på nytt **sr** Nije za ponovnu upotrebu **ru** Не пригоден для повторного использования **tr** Tekrar kullanım için değil

de Nicht erneut sterilisieren **en** Do not resterilise **fr** Ne pas restériliser **es** No esterilizar de nuevo **it** Non sterilizzare nuovamente **nl** Niet opnieuw steriliseren **da** Sterilisér ikke igen **el** Μη το αποστειρώνετε εκ νέου **pt** Não esterilizar novamente **fi** Ei saa steriloida uudelleen **sv** Sterilisera inte en gång till **et** Mitte reesteriliseerida **lv** Nesterilizēt atkārtoti **lt** Nesterilizuoti dar kartą **pl** Nie sterylizować ponownie **sk** Znovu nesterilizujte **sl** Ne ponovno sterilizirati **cs** Opakovaně nesterilizujte **hu** Nem újra sterilizálható **bg** Да не се стерилизира повторно **ro** A nu se resteriliza **hr** Ne ponovno sterilizirati **no** Ikke resteriliser **sr** Ne sterilisati ponovo **ru** не подвергать повторной стерилизации **tr** Yeniden sterilize edilemez

STERILE | EO

de Sterilisation mit Ethylenoxid **en** Sterilisation with ethylene oxide **fr** Stérilisation à l'oxyde d'éthylène **es** Esterilización con óxido de etileno **it** Sterilizzazione con ossido di etilene **nl** Sterilisatie met ethyleenoxide **da** Steriliseret med ethylenoxid **el** Αποστείρωση με οξείδιο του αιθυλενίου **pt** Esterilização com óxido de etileno **fi** Sterilointimenetelmä - etyleenioksidilla **sv** Steriliserad med etylenoxid **et** Steriliseeritud etüleenoksiidiga **lv** Sterilizēt ar etilēnoksīdu **lt** Sterilizacija etileno oksidu **pl** Wyjaławianie tlenkiem etylenu **sk** Sterilizácia s etylén oxidom **sl** Sterilizacija z etilen oksidom **cs** Sterilizace etylén oxidem **hu** Sterilizálás etilénoxiddal **bg** Стерилизация с етиленов оксид **ro** Sterilizare cu oxid de etilenă **hr** Steriliziranje sa etilenoksidom **no** Sterilisering med etylenoksid **sr** Sterilizacija sa etilenoksidom **ru** Стерилизация этиленоксидом **tr** Etilenoksit ile sterilizasyon

de Doppelples Sterilbarriersystem **en** Double sterile barrier system **fr** Système de barrière stérile double **es** Sistema de barrera estéril doble **it** Sistema a barriera sterile doppia **nl** Dubbel steriel barrièresysteem **da** Dobbelt sterilt barrieresystem **el** Διπλό σύστημα στείρου φράγματος **pt** Sistema de barreira estéril dupla **fi** Kaksinkertainen steriili sulkujärjestelmä **sv** Dubbelt sterilt barriärsystem **et** Kahekordne steriilne tõkkesüsteem **lv** Dubultā sterilās barjeras sistēma **lt** Dviguba sterili barjerinė sistema **pl** Podwójna sterylna bariera **sk** Dvojité systém sterilnej bariéry **sl** Dvojni sterilni pregradni sistem **cs** Dvojité systém sterilní bariéry **hu** Duplán sterilizált védőrendszer **bg** Двойна стерилна бариерна система **ro** Sistem de barieră sterilă dublă **hr** Dvostruki sterilni barijerni sustav **no** Dobbelt sterilbarriersystem **sr** Dupli sterilni barijerni sistem **ru** Двойная стерильная барьерная система **tr** Çift steril bariyer etkisi

de Bei beschädigter Verpackung nicht verwenden **en** Do not use if the packaging is damaged **fr** Ne pas utiliser en cas d'emballage endommagé **es** No usar cuando el empaque está dañado **it** Non usare se la confezione è danneggiata **nl** Bij beschadigde verpakking, niet gebruiken **da** Må ikke anvendes, hvis emballagen er beskadiget **el** Μη το χρησιμοποιείται σε περίπτωση που έχει υποστεί ζημιά ή συσκευασία **pt** Não utilizar se a embalagem estiver danificada **fi** Ei saa käyttää, jos pakkaus on vioittunut **sv** Använd inte vid skadad förpackning **et** Mitte kasutada, kui pakend on kahjustatud **lv** Nelietot, ja bojāts iepakojums **lt** Nenaudoti, jei pakuotė pažeista **pl** Jeżeli opakowanie jest uszkodzone – nie stosować **sk** Nepoužívať pri poškodenom balení **sl** Ne uporabljati v primeru poškodovane embalaže **cs** Je-li poškozený obal, nepoužívejte **hu** Sérült csomagolás esetén ne használja fel **bg** Да не се използва при нарушена опаковка **ro** A nu se utiliza dacă ambalajul este deteriorat **hr** Kod oštećenog omota ne koristiti **no** Ikke bruk dersom pakningen er skadet **sr** Ne upotrebljavati ako je pakovanje oštećeno **ru** Не использовать, если нарушена упаковка **tr** Hasarlı ambalajlardaki malzeme kullanılmamalıdır

de Vor direkter Sonneneinstrahlung schützen **en** Protect from direct sunlight **fr** Protéger contre les rayons directs du soleil **es** Proteger de rayos directos del sol **it** Proteggere dai raggi solari diretti **nl** Tegen direct zonlicht beschermen **da** Beskyttes mod direkte sollys **el** Προστατέψτε το από την άμεση ακτινοβολία ηλίου **pt** Proteger da incidência directa do sol **fi** Suojattava suoralta auringonpaisteelta **sv** Skydda mot direkt solstrålning **et** Kaitsta otsese päikesevalguse eest **lv** Sargāt no tiešiem saules stariem **lt** Saugoti nuo tiesioginių saulės spindulių **pl** Chronić przed bezpośrednim nasłonecznieniem **sk** Chránťe pred priamym slnečným žiarením **sl** Zaščititi pred neposrednim osončenjem **cs** Chraňte před přímými slunečnými paprsky **hu** Közvetlen napsugárzástól óvja **bg** Да се пази от пряка слънчева светлина **ro** A se feri de radiația solară directă **hr** Štititi od direktnih sunčanih zraka **no** Beskyttes mot direkte sollys **sr** Ne izlagati direktnim sunčevim zracima **ru** Защищать от прямых солнечных лучей **tr** Doğrudan gelen güneş ışığına karşı koruyunuz


de Trocken aufbewahren **en** Store in a dry place **fr** A conserver au sec **es** Mantener seco **it** Conservare in un luogo asciutto **nl** Droog bewaren **da** Opbevares tørt **el** Διατηρείστε το σε ξηρό μέρος **pt** Guardar num local seco **fi** Säilytettävä kuivassa **sv** Förvara torrt **et** Hoida kuivas **lv** Uzglabāt sausā vietā **lt** Laikyti sausai **pl** Przechowywać na suchu **sk** Uchovávať na suchom mieste **sl** Shranjevati na suhem **cs** Uchovávejte v suchu **hu** Száraz helyen tárolja **bg** Да се съхранява на сухо **ro** A se păstra la loc uscat **hr** Čuvati na suhom **no** Oppbevares tørt **sr** Čuvati na suvom **ru** Хранить в сухом месте **tr** Kuru yerde muhafaza ediniz

de Diethylhexylphthalat-(DEHP)-frei **en** Does not contain diethylhexylphthalate (DEHP)
fr sans Diéthylhexylphthalate (DEHP) **es** Libre de di(2-ethylhexil)ftalato (DEHP) **it** Senza
di-2-etilesilftalato -(DEHP) **nl** Di-ethylhexylftalaat-(DEHP)-vrij **da** Diethylhexylphtha-
lat-(DEHP)-fri **el** Ελεύθερο διεθυλεξυλ-φθαλικού εστέρα (DEHP) **pt** isento de dietilhe-
xilftalat (DEHP) **fi** Ei dietyyliheksyyliftalaattia (DEHP) **sv** Dietylhexylftalat-(DEHP)-
fritt **et** Dietüülheksüülftalaadi-(DEHP)-vaba **lv** Nesatur dietileksilflatātu (DEHP) **lt** Be
dietilheksilftalato (DEHP) **pl** Nie zawiera ftalanu dwuoktylu (DEHP) **sk** Neobsahuje di-
etylhexylftalát (DEHP) **sl** Brez dietil heksil ftalatov (DEHP) **cs** Volný diethylhexylftalát
(DEHP) **hu** Dietil-hexil-ftalát-(DEHP)-mentes **bg** Не съдържа диетилхексилфталат
(DEHP) **ro** Fără dietilhexilftalat (DEHP) **hr** Bez dietil heksilftalata-(DEHP-) **no** Die-
tylheksylftalat-(DEHP)-fri **sr** Ne sadrži dietilheksilftalat (DEHP) **ru** Не содержит
диэтилгексилфталат (DEHP) **tr** Dietil heksil-ftalat-(DEHP)-içermez

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TiLOOP® Bra

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
TiLOOP® Bra

Description

The sterile mesh implants comprise a large-pored non-absorbable prosthetic knitted fabric made using a monofilament polypropylene thread with a titanium containing coating. The mesh implants are available in various grammages with defined shapes and dimensions. The TiLOOP® breast mesh implants have been developed for use in the most diverse indications and all implantation techniques that have been applied to date in the field of reconstructive and plastic-aesthetic breast surgery. The mesh implants are designed in keeping with the requirements of the respective area of use in terms of their tensile strength and strain characteristics.

Intended use

The TiLOOP® breast mesh implants are intended for the support, reinforcement and bridging of body's own tissue.

Indication

TiLOOP® Bra is intended for the support, reinforcement and bridging of body's own tissue in reconstructive and plastic-aesthetic breast surgery.

Mode of action

The mesh structure serves as matrix for the ingrowing of stromal cells. Thereby, the adsorption of a thin fibrous tissue layer appears whose cells grow through the pores of the mesh fabric. By this effect, the mesh implant is completely enclosed by body's own tissue. Experimental investigations have shown that the implant remains soft and flexible.

The titanium containing coating of the plastic does not have an effect on radiological diagnosis procedures (X-ray, CT, MRI and sonographic examinations).

Warnings

- ☞ The packaging is to be carefully examined before use in order to ensure that it is intact. Do not use the product if the packaging should be damaged, soaked or accidentally opened or should you have doubts regarding the sterility or should the use by date have expired.
- ☞ Carefully inspect the mesh implant before use in order to ensure that it is intact. Do not use a damaged mesh implant. Damage is deemed to have been caused if the function of the mesh can no longer be warranted. Ensure that the mesh implant is not accidentally damaged during surgery. A conscious cutting of the mesh implant to size using a pair of scissors or a scalpel is permitted. When cutting the mesh implant to size, ensure that there is an adequate clearance from the seam so that it cannot be damaged. Avoid an irreversible overstretching of the mesh material. Do not use an irreversibly deformed mesh implant.
- ☞ Each product is packed separately and is supplied EO-sterilised and pyrogen-free. It is only to be used once. Do not reuse, reprocess or resterilise. A reuse, reprocessing or resterilisation of disposable products can result in a reduced performance or a loss of their functionality. A reuse of disposable products can result in exposure to pathogens such as viruses, bacteria, fungi or prions.
- ☞ It is to be taken into account in connection with all events that result in a change to the breast volume (e.g. growth phase, weight changes, pregnancies) that the mesh implant only has a restricted degree of elasticity.

Precautions

- The surgeon should be experienced with the surgical procedures and techniques. The surgeon is responsible for the implantation.
- It is to be ensured that the mesh implant is positioned without any creases if possible.

Contraindications

- Do not use in patients with an intolerance and / or allergy / allergies in connection with the product.
- Do not use in infected, contaminated or inflamed areas.

Recommended implantation procedure

General

The TiLOOP® breast mesh implant can be used in various forms of oncological and plastic-aesthetic breast surgery. The following surgical techniques are possible:

- Reconstruction after ablatio mammae (modified radical mastectomy, skin sparing mastectomy and nipple-areola-complex sparing subcutaneous mastectomy)
- Operations for tightening, elevating and lifting

Application

Step 1:

The musculus pectoralis major is dissected and divided from the thoracic wall from its caudal end towards the sternal edge. The breast implant pocket or expander pocket that is to be formed for the first time is formed cranial from medial to lateral, using the separated muscle as a covering surface for the breast implant.

Step 2:

The TiLOOP® Bra mesh implant is positioned in front of the breast implant or expander. It is attached at the free caudal and lateral end of the Musculus pectoralis major using a continuous or interrupted suture and the lower end of the mesh implant is loosely positioned behind the breast implant or expander so that it is flapped or fixed in place with a suture. It is important to cover the lower edge of the breast implant or expander with the mesh implant. When anchoring TiLOOP® Bra cranially to the muscular edge it has to be paid attention to a flexible and elastic arrangement of the mesh material.

Step 3:

The chosen breast implant or expander is inlaid below the musculus pectoralis major and placed according to the desired position.

Step 4:

The unattached part of the mesh implant that remains after cranial fixation of the suture to the edge of the Musculus pectoralis major is positioned underneath and behind the breast implant or expander in the area of the lower breast crease, where it is positioned without any creases if possible or it can be fixed in place using monofilament suturing material by means of interrupted sutures in such a way that a dislocation of the mesh implant is prevented.

Postoperative Recommendations

- The wearing of a supporting bra and a so-called Stuttgart belt is urgently recommended.
- The lifting of heavy objects, manual work and sporting activities should be avoided for 4-6 weeks.

Undesirable side-effects

- | | |
|---------------------------|---|
| • Seroma | • Infection |
| • Foreign body reaction | • Wound healing disturbance |
| • Capsular fibrosis | • Pain |
| • Post-operative bleeding | • Unsatisfactory aesthetic outcome |
| • Haematoma | • Palpable mesh structure
(if the skin is very thin) |
| • Wound dehiscence | • Temporary reactive process |
| • Necrosis | |
| • Dysaesthesia | |

Material information

Tab. 01: Overview of the materials

Mesh implant materials	Material shares that the patient comes into contact with
titanium containing coating	mainly
polypropylene	slightly (e.g. when the mesh implant is cut to size)

The titanized mesh implant does not contain diethylhexylphthalate (DEHP).

Disposal


All of the components that are supplied with the product and any mesh implants that are not used, mesh blanks and explants are to be handled and disposed of safely conform with standard hospital rules and national and regional laws and regulations. Components that are supplied together with the product could be packaging and instructions for use, for example. Special attention is to be paid to the special hazardous waste requirements in the case of infectious waste.

Miscellaneous

- The patient is to be informed of the relevant warnings, precautions, contraindications, residual risks, undesirable side-effects, and necessary follow-up.
- The patient is to be issued with the patient information leaflet (PI) and the completed implant card (IC). The PI provides the patient with general information about the mesh implant and the necessary follow-up measures. The IC informs the patient where to find further information.
- The link to the summary of safety and clinical performance (SSCP = summary of safety and clinical performance) shall be made available as soon as the European database on medical devices (Eudamed = European database on medical devices) has been completed. Until this time, the SSCP is available under the following website: www.pfmmedical.com/downloads
- All serious incidents that occur in connection with the product are to be reported to the manufacturer and the responsible national authority.
- When using the electronic version of the instructions for use, make sure that its issue date is at least as recent as the instructions for use included with the product. The current version of the electronic instructions for use can be downloaded from the manufacturer's homepage under the following link: www.pfmmedical.com/downloads

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