Titanized Mesh Implants

TiLENE® Blue Mesh

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TiLENE® Blue Mesh

Description

A large-pored non absorbable prosthetic mesh comprising monofilament polypropylene threads with an enclosing coating containing titanium and a blue stripe for orientation. The mesh exactly matches the dynamics of the body tissue as regards the tear resistance and the stretching properties. The mesh implants are available in defined forms and dimensions. These Instructions for Use are also valid for other TiLENE® Blue Mesh product variants.

Indication

TiLENE® Blue Mesh is for the surgical care or avoidance of tissue defects in all hernia indications.

TiLENE® Blue Mesh is particularly suitable for laparoscopic and open repair of direct and indirect inguinal, femoral, umbilical and incisional hernias (including pre-peritoneal and intra-peritoneal mesh positioning) and for parietal tissue strengthening and repair of the abdominal wall.

Mode of action

The purpose of the implants is to strengthen tissue or bridge defects in soft tissue. Adsorption of a thin fibrous tissue layer, whose cells grow through the pores of the mesh fabric, results in the complete enclosure of the mesh material by the patient's own body tissue. Experimental investigations have shown that the implant remains soft and flexible in this context.

Use

TiLENE® Blue Mesh is positioned as required in view of the selected and standardized surgical method. An adequate preparation is to be ensured when doing so. The size of the mesh is to be so selected that the defect is sufficiently overlapped (normally approx. 5 cm on each side). A required individual adaptation of the mesh can be carried out using a pair of scissors or a scalpel. Ensure that the mesh is positioned free from tension and creasing.

Fixation

If a fixation should be considered necessary, ensure that this is positioned approx. 1 cm from the edge of the mesh and that the distance between each of the fixation points does not exceed 2 cm. TiLENE® Blue Mesh can be fixed in place using suture material, clips, staples, clamps or adhesive. The surgeon is to select the suture material, the size of the clamps and the distances between them and the amount of adhesive with the distribution, taking the latest medical standards into account.

Warnings

- ☞ The following complications have been observed in the treatment of soft tissue defects with TiLENE® Blue Mesh mesh implants: Seroma, hematoma, recurrence, infections, visceral adhesions, mesh migration and fistula formation.
- ☞ Although the mesh material is not affected by infections, the use of non-resorbable mesh material in an infected or contaminated wound always carries the risk of fistula formation and/or mesh extrusion.
- ☞ A foreign-body reaction triggered by the titanium coating may occur in the case of a titanium allergy.
- ☞ The titanium coating of the synthetic material will not affect radiological diagnostic techniques (X-ray, CT, MRI and ultrasonic examinations).
- ☞ Do not use the product if the packaging is open or damaged or should you have doubts about the sterility or if the use by date has expired.
- ☞ Each product is packed separately and is supplied EO-sterilised and pyrogen-free. It is only to be used once. Do not reuse, reprocess or resterilise. A reuse, reprocessing or resterilisation of disposable products can result in a reduced performance or a loss of their functionality. A reuse of disposable products can result in exposure to pathogens such as viruses, bacteria, fungi or prions.
- ☞ Damaged implants must not be used. Make sure that sharp objects do not affect the mesh beyond the usual clinical application during the operation.

Precautions

- In the case of wide dissection and large-scale tissue excision, it is recommended that hemostasis be carefully checked and any remaining cavity avoided wherever possible, so as to facilitate the rapid integration of hernia strengthening.
- This product should only be used by doctors who are familiar with surgical techniques of this type. The surgeon bears the responsibility for the implantation.
- The implant should only be handled with powder-free gloves to prevent the inflammatory effect of foreign particles such as talc.
- The implants are packaged in a sterile double-wrap packaging. Carefully check that the packaging is undamaged before use. Under no circumstances should the implant be used if the packaging is damaged, open or moist. The inner packaging should only be opened directly before insertion of the implant and should exclusively be handled with sterile gloves and sterile instruments.
- Before using TiLENE® Blue Mesh mesh implants, the physician should check the size of the hernial orifice and ensure a sufficient overlap.

Contraindications

All established standard contraindications and rules for soft-tissue implants apply to the use of TiLENE® Blue Mesh mesh implants. These include the following:

- Patients in the growth period (the known extensibility of this strengthening may be insufficient for the growth of children).
- Operation on a contaminated or inflamed area.

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