Titanized Mesh Implants

TiLENE® Inguinal-L

TiLENE® Inguinal-R

de	Gebrauchsanweisung	5
en	Instructions for use	8
fr	Notice d'utilisation	11
es	Instrucciones de uso	14
pt	Instruções de utilização	17
it	Istruzioni per l'uso	20
nl	Gebruiksaanwijzing	23
sv	Bruksanvisning	26
fi	Käyttöohjeet	29
da	Brugsanvisning	32
tr	Kullanma kılavuzu	35
el	Οδηγίες χρήσεως	38
pl	Instrukcja użytkowania	41
cs	Návod k použití	44
sk	Návod na použitie	47
ru	ИНСТРУКЦИЯ ПО ИСПОЛЬЗОВАНИЮ	50
sr	Uputstvo za upotrebu	53

Erklärung der Symbole auf Etikett und Verpackung // Explanation of the symbols on label and packaging // Explication des symboles sur l'étiquette et l'emballage // Explicación de los símbolos en la etiqueta y el envoltorio // Explicação dos símbolos na etiqueta e embalagem // Illustrazione die simboli riportati su etichetta e confezione // Verklaring van de symbolen op het etiket en de verpakking // Förklaring av symboler på etiketten och förpackningen // Etiketin ja pakkauksen symbolien selitys // Forklaring af symbolerne på etiket og emballage // Etiket ve ambalaj üzerindeki sembollerin açıklamaları // Εξήγηση των συμβόλων πάνω στην ετικέτα και τη συσκευασία // Objaśnienie symboli na etykietcie i opakowaniu // Vysvětlivky k symbolům na etiketě a balení // Vysvetlivky symbolov na štítkoch a obaloch // Объяснение символов на этикетке и упаковке // Objašnjenje simbola na etiketi i na pakovanju

Hersteller // Manufacturer // Fabricant // Fabricante // Fabricante // Produttore // Fabrikant // Tillverkare // Valmistaja // Producent // Üretici // κατασκευαστής // Producent // Výrobce // Выrobca // Изготовитель // Proizvodač

Bestellnummer // Reference number // N° de commande // Número de referencia para pedidos // Número de encomenda // Numero d'ordine // Bestellnummer // Beställningsnummer // Viitenumero // Bestillingsnummer // Sipariş numarası // Αριθμός καταλόγου // Numer zamówienia // Objednací číslo // Objednávacie číslo // Номер для заказа // Broj narudbe

Chargen-Nummer // Batch number // N° de lot // Número de lote // Número do lote // Numero del lotto // Chargennummer // Satsnummer // Eränumero // Lotnummer // Yük numarası // Κωδικός παρτίδας // Kod partii // Číslo šarže // Číslo šarže // Номер партии // Serijski broj

Gebrauchsanweisung beachten // Follow Instructions for Use // Respecter les instructions du mode d'emploi // Respete las instrucciones // Cumprir as instruções de utilização // Seguire le modalità d'uso // Gebruiksaanwijzing raadplegen // Följ bruksanvisningen // Katso käyttöohje // Følg brugsanvisningen // Kullanim talimatina dikkat // Προσέξτε τις οδηγίες χρήσης // Zapoznaj się z instrukcją stosowania // Dodržujte návod k použití // Dodržujte návod na použitie // Соблюдайте инструкцию по эксплуатации // Obratiti pažnju na uputstva za upotrebu

Vor direkter Sonneneinstrahlung schützen // Protect from direct sunlight // Protéger contre les rayons directs du soleil // Proteger de rayos directos del sol // Proteger da incidência directa do sol // Proteggere dai raggi solari diretti // Tegen direct zonlicht beschermen // Skydda mot direkt solstrålning // Suojattava suoralta auringonpaisteelta // Beskyttes mod direkte sollys // Doğrudan gelen güneş ışığına karşı koruyunuz // Προστατέψτε το από την άμεση ακτινοβολία ηλίου // Chronić przed bezpośrednim nasłonecznieniem // Chraňte před přímými slunečními paprsky // Chráňte pred priamym slnečným žiarením // Защищать от прямых солнечных лучей // Ne izlagati direktnim sunčevim zracima

Trocken aufbewahren // Store in a dry place // A conserver au sec // Manter seco // Guardar num local seco // Conservare in un luogo asciutto // Droog bewaren // Förvara torrt // Säilytettävä kuivassa // Opbevares tørt // Kuru yerde muhafaza ediniz // Διατηρείστε το σε ξηρό μέρος // Przechowywać na sucho // Uchovávejte v suchu // Uchovávať na suchom mieste // Хранить в сухом месте // Čuvati na suvom

Verwendbar bis // Use by // Date limite d'utilisation // Úsese hasta // Utilizar de preferència antes de // Da usare fino a // Niet te gebruiken na // Används före // Viimeinen käytöpäivämäärä // Anvendes senest // Kullanna süresi // Ημερομηνία λήξης // Zużył przed // Použitelné do // Použitelné do // Годен до // Upotrebljivo do

Nicht zur Wiederverwendung // Do not reuse // Produit non réutilisable // No reutilizable // Não reutilizável // Non riciclabile // Niet bestemd voor hergebruik // För engångsbruk // Ei saa käyttää uudelleen // Må ikke genanvendes // Tekrar kullanım için değil // Μόνο για μία χρήση // Nie używać powtórnie // Ne k opakovanému použití // Nie k opätownému použití // Не пригоден для повторного использования // Nije za ponovnu upotrebu

Nicht erneut sterilisieren // Do not reesterilise // Ne pas restériliser // No esterilizar de nuevo // Não esterilizar novamente // Non sterilizzare nuovamente // Niet opnieuw steriliseren // Steriliserar inte en gång till // Ei saa steriloida uudelleen // Sterilisér ikke igen // Yeniden sterilize edilemez // Μη το αποστειρώνετε εκ νέου // Nie sterylizować ponownie // Opakovaně nesterilizujte // Znovu nesterilizujte // Не подвергать повторной стерилизации // Ne sterilisati ponovo

STERILE **EO**

Sterilisation mit Ethylenoxid // Sterilisation with ethylene oxide // Stérilisation à l'oxyde d'éthylène // Esterilización con óxido de etileno // Esterilização com óxido de etileno // Sterilizzazione con ossido di etilene // Sterilisatie met ethyleenoxide // Steriliserad med etylenoxid // Sterilointimenetelmä - etyleenioksidit // Steriliseret med ethylenoxid // Etilenoksit ile sterilizasyon // Αποστείρωση με οξείδιο του αιθυλενίου // Wyjaławianie tlenkiem etylenu // Sterilizace etylén oxidem // Sterilizácia s etylén oxidom // Стерилизация этиленоксидом // Sterilizacija sa etilenoksidom

Bei beschädigter Verpackung nicht verwenden // Do not use if the packaging is damaged // Ne pas utiliser en cas d'emballage endommagé // No usar cuando el empaque está dañado // Não utilizar se a embalagem estiver danificada // Non usare se la confezione è danneggiata // Bij beschadigde verpakking, niet gebruiken // Använd inte vid skadad förpackning // Ei saa käyttää, jos pakkaus on vioittunut // Må ikke anvendes, hvis emballagen er beskadiget // Hasarlı ambalajlardaki malzeme kullanılmamalıdır // Μη το χρησιμοποιείται σε περίπτωση που έχει υποστεί ζημιά ή συσκευασία // Jeżeli opakowanie jest uszkodzone - nie stosować // Je-li poškozený obal, nepoužívejte // Nepoužívajte pri poškodenom balení // Не использовать, если нарушена упаковка // Ne upotrebljavati ako je pakovanje oštećeno


Diethylhexylphthalat-(DEHP)-frei // Does not contain diethylhexylphthalate (DEHP) // sans Diethylhexylphthalate (DEHP) // Libre de di(2-ethylhexyl)ftalato (DEHP) // isento de dietilhexilftalat (DEHP) // Senza di-2-etilesilftalato -(DEHP) // Di-ethylhexylftalaat-(DEHP)-vrij // Dietylhexylftalat-(DEHP)-fritt // Ei dietyliheksyyliftalaattia (DEHP) // Diethylhexylphthalat-(DEHP)-fri // Dietil heksil-ftalat-(DEHP)-içermez // Ελεύθερο διεθυλεξυλ-φθαλικού εστέρα (DEHP) // Nie zawiera ftalanu dwuoktynu (DEHP) // Volný diethylhexylftalát (DEHP) // Neobsahuje dietylhexylftalát (DEHP) // Не содержит диэтилгексилфталат (DEHP) // Ne sadrži dietilheksilftalat (DEHP)

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
TiLENE® Inguinal-L

TiLENE® Inguinal-R

TiLENE® Inguinal-L

TiLENE® Inguinal-R

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

TiLENE® Inguinal-R

Description

Large-pore, non-absorbable prosthetic mesh made from monofilament polypropylene thread with surrounding titanized coating. The mesh is perfectly coordinated with the dynamics of body tissue in terms of resistance to tearing and elasticity. The mesh implants are available in various basis weights with defined shapes and dimensions.

Indication

Surgical mesh implants of titanized polypropylene mesh are used in the treatment of soft-tissue defects. TiLENE® Inguinal-L and TiLENE® Inguinal-R are used to repair inguinal hernias.

Mode of action

The purpose of the implants is to strengthen tissue or bridge defects in soft tissue. Adsorption of a thin fibrous tissue layer, whose cells grow through the pores of the mesh fabric, results in the complete enclosure of the mesh material by the patient's own body tissue. Experimental investigations have shown that the implant remains soft and flexible in this context.

Warnings

- ⚠ **Although the mesh material is not affected by infections, the use of non-resorbable mesh material in an infected or contaminated wound always carries the risk of fistula formation and/or mesh extrusion.**
- ⚠ **A foreign-body reaction triggered by the titanium coating may occur in the case of a titanium allergy.**
- ⚠ **The titanium coating of the synthetic material will not affect radiological diagnostic techniques (X-ray, CT, MRI and ultrasonic examinations).**
- ⚠ **Do not use the product if the packaging is open or damaged or should you have doubts about the sterility or if the use by date has expired.**
- ⚠ **Each product is packed separately and is supplied EO-sterilised and pyrogen-free. It is only to be used once. Do not reuse, reprocess or resterilise. A reuse, reprocessing or resterilisation of disposable products can result in a reduced performance or a loss of their functionality. A reuse of disposable products can result in exposure to pathogens such as viruses, bacteria, fungi or prions.**
- ⚠ **Damaged implants must not be used. Make sure that sharp objects do not affect the mesh beyond the usual clinical application during the operation.**

Precautions

- In the case of wide dissection and large-scale tissue excision, it is recommended that hemostasis be carefully checked and any remaining cavity avoided wherever possible, so as to facilitate the rapid integration of hernia strengthening.
- This product should only be used by doctors who are familiar with surgical techniques of this type. The surgeon bears the responsibility for the implantation.
- The implant should only be handled with powder-free gloves to prevent the inflammatory effect of foreign particles such as talc.
- The implants are packaged in a sterile double-wrap packaging. Carefully check that the packaging is undamaged before use. Under no circumstances should the implant be used if the packaging is damaged, open or moist. The inner packaging should only be opened directly before insertion of the implant and should exclusively be handled with sterile gloves and sterile instruments.
- Before using TiLENE® mesh implants, the physician should check the size of the hernial orifice and ensure a sufficient overlap.

Recommended Implantation procedure

TiLENÉ® Inguinal-L is used for a hernia on the left and TiLENÉ® Inguinal-R is used for a hernia on the right.

Step 01: Preparing the TiLENÉ® Inguinal-L and TiLENÉ® Inguinal-R

In order to prepare them, the edge that overlaps both of the mesh parts flush, is upturned. Using a pair of scissors, an incision is made in the larger mesh (basic mesh) below the upper smaller layer of mesh (flap) in the centre from above, approx. half of which being parallel with the suture. The end of the incision should still be concealed by the flap. When making an incision in the mesh, ensure that the clearance from the suture suffices so as to ensure that it is not damaged. The lower edge of the basic mesh can be cut so that it is conform with the anatomical situation.

Step 02: Positioning the TiLENÉ® Inguinal-L and TiLENÉ® Inguinal-R

The incision in the mesh is in a cranial position. The flap is folded over in a medial direction and the lateral part of the basic mesh is folded in a caudal direction. This part of the basic mesh that is folded over is inserted under the funiculus from medial to lateral. The lateral part of the basic mesh is folded out, followed by the flap.

The mesh should overlap medial from the edge of the rectus and the edge of the mesh should hereby extend to the symphysis. In a caudal direction, the mesh should overlap the upper pubic bone by 1-1.5 cm (especially in the area of the femoral hernia hilum) and in a cranial direction, it should generously overlap the transverse arcade. In a lateral direction, the mesh should only just extend to the Spina iliaca anterior superior.

Step 03: Fixation of the TiLENÉ® Inguinal-L and TiLENÉ® Inguinal-R

A fixation is not normally necessary. Should (in the case of major direct hernias for example), a fixation nevertheless be necessary, this should only be carried out medial at the Cooper's suspensory ligaments.

Contraindications

All established standard contraindications and rules for soft-tissue implants apply to the use of TiLENÉ® mesh implants. These include the following:



- Patients in the growth period (the known extendibility of this strengthening may be insufficient for the growth of children).
- Operation on a contaminated or inflamed area.

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