

Vascular Access Midline

primeMIDLINE®

Power Injectable, Low-Profile Midline Catheters for reliable peripheral access of selected intravenous therapies and power injection of contrast media.

www.pfmmedical.co.uk



The power Injectable **primeMIDLINE® Catheter** is designed to provide reliable peripheral access to the venous system. The **primeMIDLINE® Catheter** combines soft polyurethane material and a low-profile luer hub for patient comfort. All lumens are power injectable to eliminate the need to select the appropriate lumen and minimize interruption of work flow. The **primeMIDLINE® Catheter** offers the only high flow power injectable 3F midline catheter rated 5ml/sec.

Benefits and Features

Soft Polyurethane Material - Body softness of the material to minimise risk of damage to the inner vessel wall and provide easier insertion of the device.

High Injection Rates - Provides optimal imaging for contrast studies. 5ml/sec for 3F Single Lumen Catheter and 7ml/sec for 4F and 5F single and Dual Lumen Catheters.

Low-Profile Luer Hubs - Lays flat for patient comfort under occlusive dressings, allows for easier grasping and twisting when changing connections and resists cracking and deformation.

3F High Performance Offering - Accommodates smaller vessels and minimises thrombotic risks. Able to achieve high flow performance at 5ml/sec.

True to Size Catheter - With no Taper at the Distal End.

Pack Contents

- ▶ Midline Catheter
- ▶ Tuohy Borst Side port with Pre-Loaded Stylet
- ▶ Safety Introducer Needle 21G x 7cm with Echo Tip
- ▶ Guidewire 45cm
- ▶ Tearaway Introducer
- ▶ Syringe 10cc
- ▶ Safety Scalpel
- ▶ Tape Measure
- ▶ Neutral Needleless Valve
- ▶ Statlock® Plus Securement Device

Ordering Information

Ref -	Lumens -	French Size -	Catheter Length (cm) -	Lumen Gauge -	Max CT Rate (ml/sec) -	Max Pressure (psi) -	PU -
PFM3SML20P	Single	3F	20	19G	5	300	5
PFM4SML20P	Single	4F	20	18G	7	300	5
PFM5SML20P	Single	5F	20	17G	7	300	5
PFM5DML20P	Dual	5F	20	18G / 18G	7	300	5

Contact

Should you have any questions our Customer Solutions team will be glad to assist you.

uk@pfmmedical.com
T +44 (0) 1625 875 388
F +44 (0) 1625 875 389

Suite 3, Armcon Business Park
London Road South, Poynton
Cheshire, SK12 1LQ, England

Manufacturer
pfm medical inc. 1916 Palomar Oaks Way, Suite 150, Carlsbad, CA 92008, USA