
Macro Digital Imaging

PathLite Compact

PathLite sets the standard for capturing precise, detailed high quality images in DICOM format to assist in diagnostic reporting, publication and training.

www.pfmmedical.co.uk



PathLite Compact

A macro-imaging system designed for clinical laboratories: premium performance in a very compact and light-weight unit. The **PathLite Compact** is a purpose built system designed specifically for gross imaging within clinical laboratories. It captures images in DICOM format in line with other images captured in hospitals. It is also LIMS compatible.

Benefits

- ▶ **The small size and weight** of the Compact unit gives users the flexibility to mount it where it's needed.
- ▶ **Complete with built-in LED lighting**, PC and medical grade accessories it is designed to optimise image capture and consistency.
- ▶ **With its wide field of view**, autofocus and 20x optical zoom, this unit offers tailored performance and high resolution with maximum flexibility.
- ▶ **Ease of calibration and efficient annotation tools**, ensure accurate, reliable measurement at all zoom levels.
- ▶ **Save time and improve productivity**, it will automate and simplify tasks enabling you to focus on what you do best - Pathology.
- ▶ **Enhance accuracy and reduce the risk of errors**, the availability of high quality images removes the need for lengthy text descriptions and hand drawings, leaving less room for error.
- ▶ **Flexible and customisable**, this light-weight unit can be placed wherever you want it in the lab and at a range of heights.
- ▶ **Data is easy to access, share and is fully traceable**, Cirdan's Vivipath software platform allows images to be accessed wherever and whenever needed. They can be shared easily for teaching, research or second opinion purposes. Full traceability ensures compliance with best practice.

Details

- ▶ Comes pre-loaded with VividPath PACS Software
- ▶ Captures images in DICOM format which are exportable to JPG, LIMS compatible
- ▶ Fast-simple operation, multiple mounts available
- ▶ Weight 1.7Kg & dimensions 120mm x 140mm x 100mm
- ▶ Polarisation filters to reduce/remove glare from wet specimens
- ▶ 1-touch calibration takes under a minute, ensures reliable measurement at all zoom levels
- ▶ Selection of control interfaces : touchscreen/keyboard/mouse/foot pedal
- ▶ Full range of medical grade accessories as standard
 - 17" touch screen PC/ keyboard/ mouse
 - Barcode reader
 - 3 position foot-pedal
 - Optional micro-phone for voice recording

Ordering Information

Ref	Description
40-14	PathLite Compact

Contact

Should you have any questions our Customer Solutions team will be glad to assist you.

uk@pfmmmedical.com
T +44 (0) 1625 875 388
F +44 (0) 1625 875 389

Suite 3, Armcon Business Park
London Road South, Poynton
Cheshire, SK12 1LQ, England

Manufacturer

Cirdan, 8 Enterprise Crescent, Ballinderry Road, Lisburn, BT28 2BP, Northern Ireland