
Primera Medical Printers

Slide and Cassette Printers

www.pfmmedical.co.uk



Slide Printer

The Primera [Signature Slide Printer](#) delivers stunning print quality on printable slides or tissue cassettes in black and colour.



Benefits

The [Signature Cassette Printers](#) are designed on a proven thermal print technology that has been used for over 30 years.

The [technology has been constantly enhanced](#) over the years and as a result printers deliver stunning print quality on printable slides or tissue cassettes in black and they are the first in their class to offer colour printing.

The [Signature Printers are Windows compatible](#) and can be used with any Windows based PC or easily integrated into an LIS system.

We offer an [optional Stand-Alone solution](#) which includes a touch screen computer, medical grade keyboard as well as a 2D scanner and preinstalled PTLab Software.

We also [supply slides and cassettes](#) for the printers to ensure the best possible service to our customers.

Primera's new [Signature Slide Printer](#) can [significantly increase the efficiency of your lab](#) while helping to reduce the risk of misidentification of specimens.

It [prints directly onto slides](#), eliminating handwriting or expensive, difficult to apply xylene-resistant labels.

It is the [only colour printer on the market](#), this is of benefit as it reduced inventory and stock of different colours of consumables.

[Slide output extension bins](#) are available to enable printed slides to be collected in larger quantities

Details

- ▶ Lowest noise emission compared to all other machines on the market
- ▶ Standard Windows driver (Vista to Windows 8) with easy integration
- ▶ Completely free PTLab and PTTTemplate software with free lifetime updates
- ▶ Prints in both black and colour which is a unique feature
- ▶ Lowest energy consumption in its class 2-16 Watts (12 Volt human-safe)
- ▶ One year warranty, no questions asked
- ▶ No mandatory maintenance contracts
- ▶ Smallest footprint in its class
- ▶ Built by a printer company with 40 years' experience and more than one million machines sold in 230 countries

Ordering Information

Ref	Description	PU
78102	Signature Slide Printer - Small Footprint, LIMS Compatible, Auto Feed	1
78108	Signature Slide Printer - high pressure feature for use with other slide types.	1
78200	Medical grade 2D barcode scanner with stand for Primera PTlab TM Softwear	1
78101-1	100-Slide Output Module including 100-Slide Output container and Output Slide Ramp	1
78101-2	Additional 100-Slide Output Container	1

Cassette Printer

The Primera [Signature Cassette Printer](#) delivers stunning print quality on printable slides or tissue cassettes in black and colour.



Benefits

The Signature Cassette printer is designed on proven thermal print technology that has been used for over 30 years. The technology has been constantly enhanced over the years and as a result printers deliver stunning print quality on printable slides or tissue cassettes in black and they are the first in their class which offer colour printing. The Signature printers are Windows compatible and can be used with any Windows based PC or easily integrated into an LIS system.

We offer an optional Stand-Alone solution which includes a touch screen computer, medical grade keyboard and mouse as well as a 2D Scanner and preinstalled PTLab Software.

We also supply slides and cassettes for the printers to ensure the best possible service to our customers.

The new Signature Cassette Printer which is available now as a manual printer can be upgraded to an automated robot system with a 4x40=160 cassette magazine.

As the system is modular, the automatic magazine feeder can be added later therefore it is possible to buy the printer first and add the magazine option later.

It is the only colour printer on the market, this is of benefit as it reduced inventory and stock of different colours of consumables.

Details

- ▶ Lowest noise emission compared to all other machines on the market
- ▶ Standard Windows driver (Vista to Windows 8) with easy integration
- ▶ Completely free PTLab and PTTemplate software with free lifetime updates
- ▶ The cassette printer prints on all cassettes lid-on or lid-off.
- ▶ The printer automatically adjusts from 35 to 45 degree surface angle which is not found in competitive units
- ▶ Prints in both black and colour which is a unique feature
- ▶ Lowest energy consumption in its class 2-16 Watts (12 Volt human-safe)
- ▶ One year warranty, no questions asked
- ▶ No mandatory maintenance contracts
- ▶ Built by a printer company with 40 year's experience and more than one million machines sold in 230 countries

Ordering Information

Ref	Description	PU
78406	Signature Cassette Printer Only, with manual feed. EU Version.	1
78408	Signature Cassette Autoloader. EU Version	1
48209	PTLab PE Software (Professional Edition)	1

Label Printer

The Primera [Label Printer LX500e](#) utilizes the latest in high-resolution inkjet technology to print brilliant full-colour product labels on-demand.



Benefits

Primera's [LX500e Colour Printer](#) is a compact, full-color, label printer designed for use in labs and other medical environments. It is ideal for printing labels, name badges, wristbands and other types of narrow format output including many waterproof substrates.

[LX500e Colour Label Printers include](#) one dye-based color ink cartridge, NiceLabel Primera Edition/BarTender UltraLite design software, printer drivers for Windows 7/8/10+ and Mac ® OS X 10.6 or higher(Mac driver coming soon), starter roll of 4" x 3" High-Gloss Inkjet Label Stock, USB 2.0 cable and instructions. Requires an available USB 2.0 port.

Details

- ▶ Optional built in label cutter
- ▶ Small, compact size
- ▶ Prints up to 10cm wide
- ▶ 4800 dpi full-colour printing

Ordering Information

Ref	Description	PU
074274	LX500e Colour Label Printer 100-240 VAC, Euro Plug, CE	1
074276	LX500eC Colour Label Printer with cutter 100-240 VAC, Euro Plug, CE	1
053374	LX500e/LX500eC Colour (CMY) ink cartridge, High Yield	1

Contact

Should you have any questions our Customer Solutions team will be glad to assist you.

uk@pfmmedical.com
T +44 (0) 1625 875 388
F +44 (0) 1625 875 389

Suite 3, Armcon Business Park
London Road South, Poynton
Cheshire, SK12 1LQ, England