

---

**Primera Medical Printers**

**Slide and Cassette Printers**

---

[www.pfmmedical.co.uk](http://www.pfmmedical.co.uk)



### Slide Printer

The Primera [Signature Slide Printer](#) delivers stunning print quality on printable slides in black and colour.



### Benefits

The [Signature Slide Printers](#) are designed on a proven thermal print technology that has been used for over 30 years. The technology has been constantly enhanced over the years and as a result printers deliver stunning print quality on printable slides in black and they are the first in their class to offer colour printing.

The [Signature Printers are Windows compatible](#) and can be used with any Windows based PC or easily integrated into an LIS system.

We offer an optional [Stand-Alone solution](#) which includes a touch screen computer, medical grade keyboard as well as a 2D scanner and preinstalled PTLab Software.

We also supply slides for the printers to ensure the best possible service to our customers.

Primera's new [Signature Slide Printer](#) can significantly increase the efficiency of your lab while helping to reduce the risk of misidentification of specimens.

It prints directly onto slides, eliminating handwriting or expensive, difficult to apply xylene-resistant labels.

It is the only colour printer on the market, this is of benefit as it reduced inventory and stock of different colours of consumables.

[Slide output extension bins](#) are available to enable printed slides to be collected in larger quantities

### Details

- ▶ Lowest noise emission compared to all other machines on the market
- ▶ Standard Windows driver (Vista to Windows 8) with easy integration
- ▶ Completely free PTLab and PTemplate software with free lifetime updates
- ▶ Prints in both black and colour which is a unique feature
- ▶ Lowest energy consumption in its class 2-16 Watts (12 Volt human-safe)
- ▶ One year warranty, no questions asked
- ▶ No mandatory maintenance contracts
- ▶ Smallest footprint in its class
- ▶ Built by a printer company with 40 years' experience and more than one million machines sold in 230 countries

### Ordering Information

Ref	Description	PU
078102	Signature Slide Printer. Small Footprint, LIMS Compatible, Auto Feed.	1
078108	Signature Slide Printer. High pressure feature for use with other slide types.	1
078142	Slide Printer EVO with air puff system, display and network connector.	1
078200	Medical grade 2D barcode scanner with stand for Primera PTLab™ software.	1
078251	Spare 100-count slide holder with smartchip (slides not included)	1
SIG078101-3	Signature extension output bin (100 slides)	1

### Cassette Printer

The Primera [Signature Cassette Printer](#) delivers stunning print quality on tissue cassettes in black and colour.



### Benefits

The Signature Cassette printer is designed on proven thermal print technology that has been used for over 30 years. The technology has been constantly enhanced over the years and as a result printers deliver stunning print quality on tissue cassettes in black and they are the first in their class which offer colour printing.

The Signature printers are Windows compatible and can be used with any Windows based PC or easily integrated into an LIS system.

We offer an optional Stand-Alone solution which includes a touch screen computer, medical grade keyboard and mouse as well as a 2D Scanner and preinstalled PTLab Software.

We also supply cassettes for the printers to ensure the best possible service to our customers.

It is the only colour printer on the market, this is of benefit as it reduced inventory and stock of different colours of consumables.

### Details

---

- ▶ Lowest noise emission compared to all other machines on the market
- ▶ Standard Windows driver (Vista to Windows 8) with easy integration
- ▶ Completely free PTLab and PTTemplate software with free lifetime updates
- ▶ The cassette printer prints on all cassettes lid-on or lid-off.
- ▶ The printer automatically adjusts from 35 to 45 degree surface angle which is not found in competitive units
- ▶ Prints in both black and colour which is a unique feature
- ▶ Lowest energy consumption in its class 2-16 Watts (12 Volt human-safe)
- ▶ One year warranty, no questions asked
- ▶ No mandatory maintenance contracts
- ▶ Built by a printer company with 40 year's experience and more than one million machines sold in 230 countries

### Ordering Information

---

Ref	Description	PU
078447	Cassette Printer EVO, with manual feed only	1
078209	PTLab PE Software (Professional Edition)	1

### Contact

**Should you have any questions our Customer Services team will be glad to assist you.**

uk@pfmmmedical.com  
T +44 (0) 1625 875 388  
F +44 (0) 1625 875 389

Suite 3, Armcon Business Park  
London Road South, Poynton  
Cheshire, SK12 1LQ, England